

Nucleic acid labeling and modification

Labeled and modified nucleotides

Fluorescent Labeled nucleotides

Fluorescent nucleotides are used as building blocks for amplifications (PCR). The fluorophore should not affect the polymerase activity, while eliciting excellent brightness in downstream processes.

* NTPs are available labeled by FluoProbes® dyes at various positions using different lengths of linkers, and covering the whole UV/VIS spectrum. They are supplied ready to use in aqueous solutions, at 5 mM concentration.

FluoProbes® dyes results in nucleotides conjugates that have superior properties compared to most of other commercially available dye-labeled NTPs :

- ◆ Excellent solubility in water
- ◆ Higher signal intensity
- ◆ Better photostability
- ◆ Lower molecular weight of dye resulting in minimal steric hindrance

* All nucleotides can also be provided labeled by several other classic or commercial dyes (over 40), on a custom basis (please inquire).

FluoProbes® Labels	425A	495A	520A	532A	550A	565A	590A
λ max absorption :	436	495	525	532	554	563	594
λ max emission :	484	527	532	553	576	592	624
N6-(4-Amino)butyl-ATP	FP-FQ4920 30 µl	FP-FQ4930 30 µl	FP-FQ4940 30 µl	FP-FQ4950 30 µl	FP-FQ4960 30 µl	FP-FQ4970 30 µl	FP-FQ4980 30 µl
N6-(6-Amino)hexyl-ATP	FP-FQ5000 30 µl	FP-FQ5010 30 µl	FP-FQ5020 30 µl	FP-FQ5030 30 µl	FP-FQ5040 30 µl	FP-FQ5050 30 µl	FP-FQ5060 30 µl
8-[(4-Amino)butyl]-amino-ATP	FP-FQ4610 30 µl	FP-FQ4620 30 µl	FP-FQ4630 30 µl	FP-FQ4640 30 µl	FP-FQ4650 30 µl	FP-FQ4660 30 µl	FP-FQ4670 30 µl
8-[(6-Amino)hexyl]-amino-ATP	FP-FQ4690 30 µl	FP-FQ4700 30 µl	FP-FQ4710 30 µl	FP-FQ4720 30 µl	FP-FQ4730 30 µl	FP-FQ4740 30 µl	FP-FQ4750 30 µl
EDA-ADP	FP-FQ4810 30 µl	FP-FQ5650 20 µl	FP-FQ4820 30 µl	FP-FQ5660 20 µl	FP-FQ5670 20 µl	FP-FQ4830 30 µl	FP-FQ4840 30 µl
EDA-ATP	FP-FQ4860 30 µl	FP-FQ5740 20 µl	FP-FQ4870 30 µl	FP-FQ5750 20 µl	FP-FQ5760 20 µl	FP-FQ4880 30 µl	FP-FQ4890 30 µl
EDA-AppNHp(EDA-AMPPNP)	FP-FQ5700 20 µl	FP-FQ7270 10 µl	FP-FQ5710 20 µl	FP-FQ7280 10 µl	FP-FQ7290 10 µl	FP-FQ5720 20 µl	FP-FQ5730 30 µl
Aminoallyl-dUTP	FP-FQ5520 20 µl	FP-FQ5530 20 µl	FP-FQ5540 20 µl	FP-FQ5550 20 µl	FP-FQ5560 20 µl	FP-FQ5570 20 µl	FP-FQ5580 20 µl
5-Propargylamino-dCTP	FP-FQ5280 20 µl	FP-FQ5290 20 µl	FP-FQ5300 20 µl	FP-FQ5310 20 µl	FP-FQ5320 20 µl	FP-FQ5330 20 µl	FP-FQ5340 20 µl
FluoProbes® Labels	610A	620A	635A	647A	655A	680A	
λ max absorption :	615	619	635	645	663	680	
λ max emission :	635	643	659	669	664	700	
N6-(4-Amino)butyl-ATP	FP-FQ4990 30 µl	FP-FQ5820 20 µl	FP-FQ5830 20 µl	FP-FQ5840 20 µl	FP-FQ5850 20 µl	FP-FQ5860 20 µl	
N6-(6-Amino)hexyl-ATP	FP-FQ5070 30 µl	FP-FQ5870 20 µl	FP-FQ5880 20 µl	FP-FQ5890 20 µl	FP-FQ5900 20 µl	FP-FQ5910 20 µl	
8-[(4-Amino)butyl]-amino-ATP	FP-FQ4680 30 µl	FP-FQ5360 20 µl	FP-FQ5370 20 µl	FP-FQ5380 20 µl	FP-FQ5390 20 µl	FP-FQ5400 20 µl	
8-[(6-Amino)hexyl]-amino-ATP	FP-FQ4760 30 µl	FP-FQ5410 20 µl	FP-FQ5420 20 µl	FP-FQ5430 20 µl	FP-FQ5440 20 µl	FP-FQ5450 20 µl	
EDA-ADP	FP-FQ5680 20 µl	FP-FQ7220 10 µl	FP-FQ7230 10 µl	FP-FQ7240 10 µl	FP-FQ7250 10 µl	FP-FQ7260 10 µl	
EDA-ATP	FP-FQ5770 20 µl	FP-FQ7360 10 µl	FP-FQ7370 10 µl	FP-FQ7380 10 µl	FP-FQ7390 10 µl	FP-FQ7400 10 µl	
EDA-AppNHp(EDA-AMPPNP)	FP-FQ7300 10 µl	FP-FQ7310 10 µl	FP-FQ7320 10 µl	FP-FQ7330 10 µl	FP-FQ7340 10 µl	FP-FQ7350 10 µl	
Aminoallyl-dUTP	FP-FQ5590 20 µl	FP-FQ7160 10 µl	FP-FQ7170 10 µl	FP-FQ7180 10 µl	FP-FQ7190 10 µl	FP-FQ7200 10 µl	
5-Propargylamino-dCTP	FP-FQ5350 20 µl	FP-FQ7110 10 µl	FP-FQ7120 10 µl	FP-FQ7130 10 µl	FP-FQ7140 10 µl	FP-FQ7150 10 µl	

Nucleic acid labeling and modification

Labeled and modified nucleotides

Label-nucleotide spacers

The number "11" is the number of atoms in the linker between biotin and dUTP. The linker length affects the incorporation efficiency of the biotin-dUTP probe into DNA using DNA polymerases, and it also affects biotin/avidin or biotin/streptavidin interaction. In general, the shorter the linker, the more efficiently the biotin-dUTP is incorporated into DNA by DNA polymerases. On the other hand, the longer the linker, the better biotin can interact with avidin or streptavidin.

Biotinylated nucleotides

Biotin is a versatile label, allowing color choice and high sensitivity through the use of various signal amplification system. Biotinylated nucleotides can be incorporated into nucleic materials by Polymerase Chain Reaction (PCR). dUTP serves for DNA amplifications. UTP serves for RNA amplifications. The resulting biotin-containing DNA or RNA can be subsequently labeled with a fluorescent labeled (strept)avidins or detected with anti-biotin antibodies.

Biotin-dUTP can be enzymatically incorporated into DNA via nick translation, random priming, or 3'-end terminal labeling.

Biotin-11-dUTP

(Biotin-11-2'-deoxyuridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 886.4

Description	Cat.#	Qty
Biotin-11-dUTP	FP-AM539A	50 µl
Biotin-11-dUTP Lyophilized powder	FP-AY489A	50 µg

Biotin-16-dUTP

(Biotin-16-2'-deoxyuridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 971.5

Description	Cat.#	Qty
Biotin-16-dUTP	FP-AM541A	50 µl
Biotin-16-dUTP, Lyophilized powder	FP-AM542A	50 µg

Biotin-20-dUTP

(Biotin-20-2'-deoxyuridine-5'-triphosphate, tetralithium salt)

MW : 1020.54

Description	Cat.#	Qty
Biotin-20-dUTP	FP-AM543A	50 µg

Biotin-11-UTP

(Biotin-11-uridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 902.5

Description	Cat.#	Qty
Biotin-11-UTP	FP-AM550A	50 µl

Biotin-16-UTP

(Biotin-16-uridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 987.52

Description	Cat.#	Qty
Biotin-16-UTP	FP-AM551A	50 µl

Biotin-EDA-ATP

2'/3'-O-(2-(Biotinyl-butyl-carbamoyl)-ethyl-carbamoyl)-Adenosine-5'-triphosphate, Na salt

MW : 815.53 (Anion)

Description	Cat.#	Qty
Biotin-EDA-ATP	FQ5630	20 u
	FQ5331	10 u

Biotin-EDA-AppNHp (Biotin-EDA-AMPPNP)

2'/3'-O-(2-(Biotinyl-butyl-carbamoyl)-ethyl-carbamoyl)-Adenosine-5'-[β,γ -imido]triphosphate, Na salt

MW : 814.55 (Anion)

Description	Cat.#	Qty
Biotin-EDA-AppNHp (Biotin-EDA-AMPPNP)	FQ2831	10 u
	FQ2830	50 u

Related Products :
Aminated nucleotides page D.137

Amine labeled nucleotides

Amine labeled nucleotides are used to link nucleotides to amine-reactive probes or to stationary phases for affinity chromatography.

N⁶-(4-Amino)butyl-ATP

N⁶-(4-Amino)butyl-adenosine-5'-triphosphate, Sodium salt
MW : 575.28 (Anion)

Description	Cat.#	Qty
N ⁶ -(4-Amino)butyl-ATP	FQ2960	50 u
	FQ2961	250 u

N⁶-(6-Amino)hexyl-ATP

N⁶-(6-Amino)hexyl-adenosine-5'-triphosphate, Sodium salt
MW : 603.33 (Anion)

Description	Cat.#	Qty
N ⁶ -(6-Amino)hexyl-ATP	FQ2970	50 u
	FQ2971	250 u

Ito, *et al.*, "Modified nucleic acid for systematic evolution of RNA ligands by exponential enrichment.", *J. Bioact. Compat. Pol.*, 13 (2), 114 (1998)

Trayer, "Affinity chromatography of some adenosine phosphate-requiring systems", *Biochem. Soc. T.*, 2 (6), 1302 (1974)

8-[(4-Amino)butyl]-amino-ATP

8-[(4-Amino)butyl]-amino-adenosine-5'-triphosphate, Sodium salt
MW : 590.32 (Anion)

Description	Cat.#	Qty
8-[(4-Amino)butyl]-amino-ATP	FQ2630	50 u
	FQ2631	250 u

8-[(6-Amino)hexyl]-amino-ATP

8-[(6-Amino)hexyl]-amino-adenosine-5'-triphosphate, Sodium salt
MW : 618.35 (Anion)

Description	Cat.#	Qty
8-[(6-Amino)hexyl]-amino-ATP	FQ2640	50 u
	FQ2641	250 u

EDA-ADP

2'/3'-O-(2-Aminoethyl-carbamoyl)-Adenosine-5'-diphosphate, Sodium salt
MW: 511.28 (Anion)

Description	Cat.#	Qty
EDA-ADP	FQ4800	30 u
	FQ4801	150 u

Oiwa, *et al.*, "Comparative Single-Molecule and Ensemble Myosin Enzymology: Sulfoindocyanine ATP and ADP Derivatives.", *Biophys. J.* 78:3048 (2000)

Jameson, *et al.*, "Fluorescent analogs: Synthesis and Applications.", *Methods in Enzymology*, 278, 363 (1997)

EDA-ATP

2'/3'-O-(2-Aminoethyl-carbamoyl)-Adenosine-5'-triphosphate, Sodium salt
MW : 590.25 (Anion)

Description	Cat.#	Qty
EDA-ATP	FQ4850	30 u
	FQ4851	150 u

Oiwa, *et al.*, "Comparative Single-Molecule and Ensemble Myosin Enzymology: Sulfoindocyanine ATP and ADP Derivatives.", *Biophys. J.* 78:3048 (2000)

Jameson, *et al.*, "Fluorescent analogs: Synthesis and Applications.", *Methods in Enzymology*, 278, 363 (1997)

Nucleic acid labeling and modification

Labeled and modified nucleotides

γ -Aminophenyl-ATP

Adenosine-5'-[γ -(4-aminophenyl)]triphosphate, Sodium salt

MW : 595.27 (Anion)

Description	Cat.#	Qty
γ -Aminophenyl-ATP	FQ5120	3 mg
	FQ5121	15 mg

Haystead, *et al.*, "Gamma-phosphate-linked ATP-Sepharose for the affinity purification of protein-kinases -rapid purification to homogeneity of skeletal-muscle mitogen-activated protein-kinase kinase.", *Eur. J. Biochem.*, 214 (2), 459 (1993)

Trayer, *et al.*, "Preparation of adenosine nucleotide derivatives suitable for affinity chromatography.", *Biochem. J.* 139 (3):609 (1974)

5-Propargylamino-Dctp

5-Propargylamino-2'-deoxy-cytidine-5'-triphosphate, Sodium salt

MW : 517.19 (Anion)

Description	Cat.#	Qty
5-Propargylamino-Dctp	FQ4580	30 μ
	FQ4581	150 μ

Hobbs, "Palladium-catalyzed synthesis of alkynylamino nucleosides - a universal linker for nucleic-acids.", *J. Org.Chem.*, 54 (14), 3420 (1989)

Related products :

Description	Cat.#	Qty
EDA-AppNHp (EDA-AMPPNP)	FQ5690	20 μ

Aminoallyl nucleotides

Aminoallyl-nucleotides can be enzymatically incorporated into nucleic materials, by Polymerase Chain Reaction (PCR). dUTP serves for DNA amplifications. The resulting amine-containing DNA or RNA can be subsequently labeled with a fluorescent dye, biotin or other haptens via conventional peptide coupling methods. This two-step method for labeling nucleic acids is considerably more economical than the one-step method using a prelabeled nucleotide.

5-Aminoallyl-DUTP

Sodium salt

(5-(3-aminoallyl)-2'-deoxyuridine-5'- triphosphate, trisodium salt, AA-dUTP)

MW : 589.17

Description	Cat.#	Qty
5-Aminoallyl-dUTP	FP-AK218A	100 μ l (10 mM in TE buffer)
	FP-AY490A	1 mg (Powder)

5-Aminoallyl-UTP

Sodium salt

(5-(3-aminoallyl)uridine-5'-triphosphate, trisodium salt, AA-UTP)

MW : 605.17

Description	Cat.#	Qty
5-Aminoallyl-UTP	FP-AM555A	100 μ l (10 mM in TE buffer)
	FP-AM558A	1 mg (Powder)

D.138 Proc Natl Acad Sci U S A 90, 4206(1993)
J. Clin Microbiol 29, 583(1991)

Nucleic acid labeling and modification

Labeled and modified nucleotides

Mant-labeled oligonucleotides

Mant fluorophore (N-Methyl-anthraniloyl) is compact so, when attached on ribose ring it causes minimal perturbation of nucleotide-protein interactions. It is however sensitive to environment, notably to conformational changes of nucleotide-binding proteins. Mant-labeled nucleotides offer valuable tools to investigate the structure, protein-ligand interactions, and enzymatic activity of nucleotide-binding proteins (as myosin with ATP, P21 with GTP).

Literature for MAnt-ATP (also available for other Mant nucleotides).

-Cheng JQ, Jiang W, Hackney DD 1998; Interaction of mant-adenosine nucleotides and magnesium with kinesin ; *Biochemistry* 37, 5289-5295

-Churchich 1997 ; Conformational changes at the nucleotide binding of GroEL induced by binding of protein substrates – Luminescence studies ; *J.Biol. Chem.*


272(32) : 19645

-Malnasi-Csizmadia A, Woolley RJ, Bagshaw CR. 2000; Resolution of conformational states of dictyostelium myosin II motor domain using tryptophan (W501)

mutant s: implications for the open-closed transition identified by crystallography. ; *Biochemistry* 39, 16135-16146

-Shaffer J, Sun G, Adams JA. 2001 ; Nucleotide release and associated conformational changes regulate function in the COOH-terminal Src kinase, *Csk*; *Biochemistry*

40, 11149-11155



NPE-caged-ATP releases ATP generally on a millisecond time scale upon flash photolysis with near-UV light.

N⁶-[4-(Mant-Amino)]butyl-ATP

N⁶-[4-((N-methyl-anthraniloyl)-amino)]butyl-adenosine-5'-triphosphate, Sodium salt

MW : 708.43 (Anion)

$\lambda_{exc} / \lambda_{em}$: 335 / 440 nm

λ_{max} : 280/335 nm ; ϵ : 18 000/2 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
N ⁶ -[4-(Mant-Amino)]butyl-ATP	FP-FQ6611	20 u
	FP-FQ6610	100 u

N⁶-[6-(Mant-Amino)]hexyl-ATP

N⁶-[6-((N-methyl-anthraniloyl)-amino)]hexyl-adenosine-5'-triphosphate, Sodium salt

MW : 736.48 (Anion)

$\lambda_{exc} / \lambda_{em}$: 335 / 440 nm

λ_{max} : 280/335 nm ; ϵ : 18 000/2 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
N ⁶ -[6-(Mant-Amino)]hexyl-ATP	FP-FQ6621	20 u
	FP-FQ6620	100 u

8-[4-(Mant-Amino)]butyl-ATP

8-[(4-(N-methyl-anthraniloyl)-amino)butyl]-amino-adenosine-5'-triphosphate, Sodium salt

MW : 723.44 (Anion)

$\lambda_{exc} / \lambda_{em}$: 335 / 440 nm

λ_{max} : 280/335 nm ; ϵ : 18 000/2 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
8-[4-(Mant-Amino)]butyl-ATP	FP-FQ6461	20 u
(MABA-ATP)	FP-FQ6460	100 u

Groemping, et al., "Regulation of ATPase and chaperone cycle of DnaK from *Thermus thermophilus* by the nucleotide exchange factor GrpE.", *J. Mol. Biol.*, 305 (5):1173(2001)

Weikl, et al., "C-terminal regions of Hsp90 are important for trapping the nucleotide during the ATPase cycle.", *J. Mol Biol.* 303 (4):583(2000)

8-[6-(Mant-Amino)]hexyl-ATP

8-[(6-(N-methyl-anthraniloyl)-amino)hexyl]-amino-adenosine-5'-triphosphate, Sodium salt

MW : 751.50 (Anion)

$\lambda_{exc} / \lambda_{em}$: 335 / 440 nm

λ_{max} : 280/335 nm ; ϵ : 18 000/2 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
8-[6-(Mant-Amino)]hexyl-ATP	FP-FQ6471	20 u
(MAHA-ATP)	FP-FQ6470	100 u

1 u = 1 μ l at 10 mM

Nucleic acid labeling and modification

Labeled and modified nucleotides

Mant-EDA-ATP

2'/3'-[(2-(N-methyl-anthraniloyl)-amino)ethyl-carbamoyl]-adenosine-5'-triphosphate, Sodium salt

MW : 723.40 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 255/355 nm ; ϵ : 23 300/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-EDA-ATP	FP-FQ5811	20 u
	FP-FQ5810	100 u

Mant-ADP

2'/3'-O-(N-Methyl-anthraniloyl)-adenosine-5'-diphosphate, Triethylammonium salt

MW : 558.32 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 255/355 nm ; ϵ : 23 300/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-ADP	FP-FQ2061	150 u
	FP-FQ2060	750 u

Bujalowski et al, "Kinetic mechanism of nucleotide cofactor binding to Escherichia coli replicative helicase DnaB protein. stopped-flow kinetic studies using fluorescent, ribose-, and base-modified nucleotide analogues.", *Biochemistry* 39:2106(2000)
Cheng et al, "Interaction of mant-adenosine nucleotides and magnesium with kinesin.", *Biochemistry* 37:5288 (1998)

Mant-ATP

2'/3'-(N-Methyl-anthraniloyl)-adenosine-5'-triphosphate, Triethylammonium salt

MW : 637.30 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 255/355 nm ; ϵ : 23 300/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-ATP	FP-FQ2071	150 u

Tung-Chung Mou, et al, "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(NMethylanthraniloyl) guanosine 5'-Triphosphate.", *J. Biol. Chem.* 280 (8):7253(2005)
Booth, et al, "Analysis of the properties of the N-terminal nucleotide-binding domain of human P-glycoprotein.", *Biochemistry* 39:5518(2000)
Thoenges, et al., "Tight binding of bulky fluorescent derivatives of adenosine to the low affinity E2ATP site leads to inhibition of Na⁺/K⁺-ATPase. Analysis of structural requirements of fluorescent ATP derivatives with a Koshland-Nemethy-Filmer model of two interacting ATP sites.", *J. Biol. Chem.* 274:1971(1999)

Mant-dATP

3'-(N-Methyl-anthraniloyl)-2'-deoxy-adenosine-5'-triphosphate, Triethylammonium salt

MW : 621.31 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 255/355 nm ; ϵ : 23 300/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-dATP	FP-FQ2901	50 u
	FP-FQ2900	250 u

Tung-Chung Mou, et al, "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(NMethylanthraniloyl)- guanosine 5'-Triphosphate.", *J. Biol. Chem.* 280 (8):7253(2005)
Churchich, "Binding of a fluorescent nucleotide analog to hsc70 - the effect of peptide protein interactions on the luminescence properties of the probe.", *Eur. J. Biochem.* 231 (3):736(1995)

Mant-GDP

2'/3'-O-(N-Methyl-anthraniloyl)-guanosine-5'-diphosphate, Triethylammonium salt

MW : 574.32 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 252/355 nm ; ϵ : 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-GDP	FP-FQ2081	150 u
	FP-FQ2080	750 u

Thanbichler, et al. "Kinetics of the interaction of translation factor SelB from Escherichia coli with guanosine nucleotides and selenocysteine insertion sequence RNA.", *J. Biol. Chem.* 275:20458(2000)
Murthy, et al, "Nucleotide binding by the erythrocyte transglutaminase/Gh protein, probed with fluorescent analogs of GTP and GDP.", *Proc. Natl. Acad. Sci. USA* 97:7744 (2000)

1 u = 1 μ l at 10 mM

Mant-dGDP

3'-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-diphosphate, Triethylammonium salt

MW : 558.33 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 252/355 nm ; ϵ : 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-dGDP	FP-FQ4901	30 u
	FP-FQ4900	150 u

Nomanbhoy, et al., "Investigation of the GTP-binding GTPase cycle of Cdc42Hs using extrinsic reporter group fluorescence.", *Biochemistry* 35 (14):4602 (1996).

Leonard, et al., "Investigation of the GTP-binding GTPase cycle of cdc42hs using fluorescence spectroscopy.", *Biochemistry* 33 (40):12323. (1994)

Mant-GTP

2'/3'-O-(N-methyl-anthraniloyl)-guanosine-5'-triphosphate, Triethylammonium salt

MW : 653.30 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 255/355 nm ; ϵ : 23 300/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-GTP	FP-FQ2091	150 u
	FP-FQ2090	750 u

Tung-Chung Mou, et al, "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(NMethylanthraniloyl)-guanosine 5'-Triphosphate.", *J. Biol. Chem.* 280 (8):7253 (2005) : the mechanism of discrimination between guanosine and adenosine nucleotides.", *Biochemistry* 34:593(1995)

Mant-dGTP

3'-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-triphosphate, Triethylammonium salt

MW : 637.30

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 252/355 nm ; ϵ : 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-dGTP	FP-FQ2911	10 u
	FP-FQ2910	50 u

Mant-GppNHp

2'/3'-O-(N-Methyl-anthraniloyl)-guanosine-5'-[(β , γ)-imido]triphosphate, Triethylammonium salt

MW : 651.31 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 252/355 nm ; ϵ : 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-GppNHp (MAnt-GMPPNP)	FP-FQ2551	10 u
	FP-FQ2550	50 u

Diebold, et al. "Molecular basis for Rac2 regulation of phagocyte NADPH oxidase", *Nature Immunol.* 2:211(2001)

Graham, et al., "The conserved arginine in rho-GTPase-activating protein is essential for efficient catalysis but not for complex formation with Rho.GDP and aluminum fluoride.", *Biochemistry* 38:985(1999)

Mant-dGppNHp

3'-O-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-[(γ , β)-imido]triphosphate, Triethylammonium salt

MW : 636.32 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} 252/355 nm ; ϵ 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-dGppNHp (MAnt-dGMPPNP)	FP-FQ2541	10 u
	FP-FQ2540	50 u

Scheidig, et al., "X-ray crystal-structure analysis of the catalytic domain of the oncogene product p21(H-Ras) complexed with caged GTP and mant dGppNHp.", *J. Mol. Biol.* 253 (1):132.(1995)

Nucleic acid labeling and modification

Labeled and modified nucleotides

Mant-GTPgammaS

2'/3'-O-(N-Methyl-anthraniloyl)-guanosine-5'-(γ -thio)-triphosphate, Triethylammonium salt

MW : 669.36 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} 252/355 nm ; ϵ 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-GTPgammaS	FP-FQ2921	10 u
	FP-FQ2920	50 u

Remmers, "Detection and quantitation of heterotrimeric G proteins by fluorescence resonance energy transfer.", *Anal. Biochem.* 257 (1):89(1998)

Remmers, et al. "Interdomain interactions regulate GDP release from heterotrimeric G proteins.", *Biochemistry* 38(42):13795(1999)

Mant-XDP

2'/3'-O-(N-Methyl-anthraniloyl)-xanthosine-5'-diphosphate, Triethylammonium salt

MW : 575.32 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} 254/355 nm ; ϵ 19 000/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-XDP	FP-FQ2941	10 u
	FP-FQ2940	50 u

Mant-XppNHp

2'/3'-O-(N-Methyl-anthraniloyl)-xanthosine-5'-[(γ , β)-imido]triphosphate, Triethylammonium salt

MW : 653.30 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} 254/355 nm ; ϵ 19 000/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-XppNHp (mant-XMPPNP)	FP-FQ5181	5 u
	FP-FQ5180	25 u

Rodnina, et al., "Codon-dependent conformational change of elongation-factor tu preceding GTP hydrolysis on the Ribosome", *EMBO J.* 14 (11):2613(1995)

Hazlett, et al., "Solution dynamics of p21(Ras) proteins bound with fluorescent nucleotides - a time-resolved fluorescence study.", *Biochemistry-US* 32 (49):13575.(1993)

Mant-XTP

2'/3'-O-(N-Methyl-anthraniloyl)-xanthosine-5'-triphosphate, Triethylammonium salt

MW : 654.28 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} 254/355 nm ; ϵ 19 000/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-XTP	FP-FQ2951	10 u
	FP-FQ2950	50 u

Tung-Chung Mou, et al., "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(N-Methylanthraniloyl)- guanosine 5'-Triphosphate.", *J. Biol. Chem.* 280 (8):7253(2005)

Mant-AppNHp

2'/3'-O-(N-Methyl-anthraniloyl)-Adenosine-5'-[(β , γ)-imido]triphosphate, Triethylammonium salt

MW : 636.32 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} 255/355 nm ; ϵ 23 300/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-AppNHp (MAnt-AMPPNP)	FP-FQ2531	10 u
	FP-FQ2530	50 u

Bujalowski, et al., "Kinetic mechanism of nucleotide cofactor binding to Escherichia coli replicative helicase DnaB protein. Stopped-flow kinetic studies using fluorescent, ribose-, and base-modified nucleotide analogues.", *Biochemistry* 39:2106(2000)

Moore, et al. "Kinetic mechanism of adenine nucleotide binding to and hydrolysis by the Escherichia coli Rep monomer. 1. Use of fluorescent nucleotide analogues.", *Biochemistry* 33:14550(1994)

Nucleic acid labeling and modification

Labeled and modified nucleotides

Mant-ITP γ S

2'/3'-O-(N-Methyl-anthraniloyl)-inosine-5'-(γ -thio)-triphosphate, Sodium salt

MW : 654.35 (Anion)

Description	Cat.#	Qty
Mant-ITP γ S	FP-FQ2931	10 u
	FP-FQ2930	50 u

Mant-N6-Methyl-ATP

2'/3'-O-(N-Methyl-anthraniloyl)-N6-methyl-adenosine-5'-triphosphate, Sodium salt

MW : 650.32 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 255/355 nm ; ϵ : 23 300/5 800 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
Mant-N6-Methyl-ATP	FP-FQ4910	30 u
	FP-FQ4911	150 u

NPE-caged-Mant-dGTP

3'-O-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-triphosphate, P3-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt

MW : 786.45 (Anion)

$\lambda_{exc} / \lambda_{em}$: 355 / 448 nm

λ_{max} : 252/355 nm ; ϵ : 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
NPE-caged-Mant-dGTP	FP-FQ3021	10 u
	FP-FQ3020	50 u

Nucleic acid labeling and modification

Labeled and modified nucleotides

Technical tip


"Caged compounds" serve as photo-induced precursors of biological effector compounds, i.e. nucleotides and Ca^{2+} chelators. Combined with reactive fluorescent probes they have been used to label biological molecules ; by FT-IR coupled with flash photolysis and by fast time-resolved UV-Visible spectroscopy coupled with laser flash photolysis.

Applications :

- ◆ Kinetic studies of release and export exchanges of nucleotides (mitochondria)
- ◆ Phosphorylation activation of enzymes (Ca^{2+} ATPase)

1) Electrical currents associated with nucleotide transport by the reconstituted mitochondrial ADP/ATP carrier ; Nickolay B. et al ; Proc. Natl. Acad. Sci. USA, Vol. 93, pp. 664-668 (1996) Article.

2) ATP-Induced phosphorylation of the sarcoplasmic reticulum Ca^{2+} ATPase : Molecular interpretation of infrared difference spectra ; A. Barth, et al (1998) ; Biophys. J. 75 538-544 / Abstract.



Capacitive currents generated by the reconstituted AAC under UV flashes with loaded (A and C) and unloaded (B and D) vesicles in the presence of caged ATP (c-ATP) or caged ADP (c-ADP), and inhibitor of CAT and BKA.

Caged nucleotides

Caged nucleotides contain a photolabile group that allows monitoring of fast kinetics especially in time-resolved X-ray crystallography : Flash photolysis of the caging group by ultraviolet light leads to a rapid and highly localized release of the biologically active nucleotide from the inactive caged compound.

Caged nucleotides are increasingly useful for :

- ◆ Identification of drug-targets and its binding site
- ◆ Determination of the affinity and selectivity of the drug-target interaction.
- ◆ Determination of the conformation of the binding pocket in structure-based drug design.

NPE-caged-ATP

Adenosine-5'-triphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt
MW : 653.30 (Anion)

Description	Cat.#	Qty
NPE-caged-ATP	FP-FQ2121	150 u
	FP-FQ2120	750 u

Scheirlinckx, et al., "Monitoring of secondary and tertiary structure changes in the gastric H⁺/K⁺-ATPase by infrared spectroscopy", Eur. J. Biochem. 268 (13):3644(2001)

Barth, et al., "Substrate binding and enzyme function investigated by infrared spectroscopy", FEBS Lett. 477:151(2000)

DMB-caged-ATP

Adenosine-5'-triphosphate, P³-(1-(3',5'-dimethoxyphenyl)-2-oxo-2-phenyl-ethyl)-ester, Triethylammonium salt
MW : 758.44 (Anion)

Description	Cat.#	Qty
DMB-caged-ATP	FP-FQ3001	10 u
	FP-FQ3000	50 u

Sokolov, et al., «Fast transient currents in Na,K-ATPase induced by ATP concentration jumps from the P-3-[1-(3',5'-dimethoxyphenyl)-2-phenyl-2-oxoethyl ester of ATP.», Biophys. J. 74 (5):2285(1998)

Sokolov, et al., «Fast transient currents in the Na,K-ATPase induced by ATP concentration jump experiments from DMB caged ATP.», Biophys. J. 72 (2):242.

NPE-caged-GTP

Guanosine-5'-triphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt
MW : 669.30

Description	Cat.#	Qty
NPE-caged-GTP	FP-FQ2421	100 u
	FP-FQ2420	500 u

Allin, et al., "Ras catalyzes CTP hydrolysis by shifting negative charges from gamma- to beta-phosphate as revealed by time-resolved FTIR difference spectroscopy.", Biochemistry-US 40 (10):3037(2001)

Scheidig, et al., "The pre-hydrolysis state of p21(ras) in complex with GTP : new insights into the role of water molecules in the GTP hydrolysis reaction of ras-like proteins.", Struct. Fold Des. 7 (11):1311(1999)

NPE-caged-mant-dGTP

3'-O-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-triphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt
MW : 786.45 (Anion)

$\lambda_{\text{exc}} / \lambda_{\text{em}}$: 355 / 448 nm

λ_{max} : 252/355 nm ; ϵ : 22 600/5 700 mol⁻¹ cm⁻¹

Description	Cat.#	Qty
NPE-caged-mant-dGTP	FP-FQ3021	10 u
	FP-FQ3020	50 u

Nucleic acid labeling and modification

Labeled and modified nucleotides

NPE-caged-AppNHp

Adenosine-5'-[(β,γ)-imido]triphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, triethylammonium salt
MW : 652.32 (Anion)

Description	Cat.#	Qty
NPE-caged-AppNHp (NPE-caged-AMPPNP)	FP-FQ2411	100 u
	FP-FQ2410	500 u

NPE-caged-GpCpp

Guanosine-5'-[(β,γ)-methylene]triphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt
MW : 667.33 (Anion)

Description	Cat.#	Qty
NPE-caged-GpCpp (NPE-caged-GMPCPP)	FP-FQ3011	10 u
	FP-FQ3010	50 u

NPE-caged-XDP

Xanthosine-5'-diphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt
MW : 591.32 (Anion)

Description	Cat.#	Qty
NPE-caged-XDP	FP-FQ3031	10 u
	FP-FQ3030	50 u

NPE-caged-XppNHp

Xanthosine-5'-[(β,γ)-imido]triphosphate, P³-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt
MW : 669.30 (Anion)

Description	Cat.#	Qty
NPE-caged-XppNHp (NPE-caged-XMPPNP)	FP-FQ5201	5 u
	FP-FQ5200	25 u

Nucleic acid labeling and modification

Labeled and modified nucleotides

BrDU nucleotides

BrDU is an analog of thymidine, used alternatively to radioactive thymidine to track DNA synthesis by (sub)populations of cells. It is added to the culture at high concentration. Newly synthesized DNA thymidine is (partially) replaced by BrdU. After fixation and permeabilization of the cells the BrdU can be shown using a fluorochrome labeled antibody against the BrdU. Using DNase during the staining period the BrdU becomes more available to the antibody without influencing the fluorescence or the structure of the cell.

BrdU

5-Bromo-2'-deoxyuridine

MW : 307.1

BrdU can be incorporated into DNA during cell division and subsequently detected by a anti-BrdU antibody. The probe can be used to study cell-cycle kinetics.

Description	Cat.#	Qty
BrdU (5-Bromo-2'-deoxyuridine)	FP-18538A	100 mg

Methods Cell Biol 41, 297(1994).

BrUTP

5-Bromouridine-5'-triphosphate

10 mM in TE buffer

MW : 628.99

BrUTP can be enzymatically incorporated into RNA.

Description	Cat.#	Qty
BrUTP	FP-T8200A	25 µl

Literature : Nakayama C, et al., "Utilizations of various uridine 5'-triphosphate analogues by DNA-dependent RNA polymerases I and II purified from liver nuclei of the cherry salmon (*Oncorhynchus masou*)", J Biochem (Tokyo) 96, 1501(1984)

Br-dUTP

5-Bromo-2'-deoxyuridine-5'-triphosphate,

10 mM in TE buffer

5-Bromo-dUTP is widely used in TUNEL assay to detect apoptosis cells(1,2) and it is also a good substrate for reverse transcriptase.

Description	Cat.#	Qty
Br-dUTP	FP-T8199A	25 µl

Li X., et al., "Detection of apoptosis and DNA replication by differential labeling of DNA strand breaks with fluorochromes of different color.", Exp Cell Res, 222, 28(1996). Li X., et al., "Labelling DNA strand breaks with BrdUTP. Detection of apoptosis and cell proliferation.", Cell Prolif, 28, 571 (1995). Rytting AS, et al., "Colorimetric capture assay for human-immunodeficiency-virus-I reverse transcriptase activity", Biotechnol Appl Biochem, 29, 241(1999)

Ekstrand DH, et al., "A sensitive assay for the quantification of reverse transcriptase activity based on the use of carrier-bound template and non-radioactive-product detection, with special reference to human-immunodeficiency-virus isolation", Biotechnol Appl Biochem, 23, 95(1996)

5Br-dUDP

5-Bromo-2'-deoxy-uridine-5'-diphosphate, Sodium salt

MW : 465.04 (Anion)

10 mM solution

Description	Cat.#	Qty
5Br-dUDP	FQ4570	30 µl
	FQ4571	150 µl

Cyclic nucleotides

Guanosine-3',5'-cyclic monophosphate

Guanosine-3',5'-cyclic monophosphate, Sodium salt

MW : 343.19 (Anion)

Description	Cat.#	Qty
cGMP	FQ2840	50 mg
	FQ2841	250 mg

Nisoli, et al. «Mitochondrial biogenesis in mammals: the role of endogenous nitric oxide.», Science, 299 (5608), 896 (2002). Rybalkin, et al, «PDE5 is converted to an activated state upon cGMP binding to the GAF A domain.», EMBO J., 22 (3), 469 (2003). Coates, et al., «Antagonistic pathways in neurons exposed to body fluid regulate social feeding in *Caenorhabditis elegans*.», Nature, 419(6910), 925 (2002)

Related products :

Description	Cat.#
PCR Cleanup Column Kit	T66380
Proteinase K	718960
TrueBlue PCR cloning kit	R58140

Nucleoside bisphosphates

Adenosine-3',5'-bisphosphate

Adenosine-3',5'-bisphosphate, Triethylammonium salt

MW : 425.18 (Anion)

Description	Cat.#	Qty
Adenosine-3',5'-bisphosphate (pAp)	FQ2281	100 u

Jacobson, et al., "Ribose modified nucleosides and nucleotides as ligands for purine receptors.", *Nucleos. Nucleot. Nucl.*, 20 (4-7), 333 (2001)

Jung, et al., "Platelet collagen receptor integrin alpha(2)beta(1) activation involves differential participation of ADP-receptor subtypes P2Y1 and P2Y12 but not intracellular calcium change.", *Eur. J. Biochem.*, 268 (12), 3513 (2001)

Sheng, et al., "Bacterial expression, purification, and characterization of rat hydroxysteroid sulfotransferase STa.", *Protein Express. Purif.*, 21 (1), 235 (2001)

Adenosine-2',5'-bisphosphate

Adenosine-2',5'-bisphosphate, Triethylammonium salt

MW : 425.18 (Anion)

Description	Cat.#	Qty
Adenosine-2',5'-bisphosphate	FQ2271	100 u
	FQ2270	500 u

Toth-Zsomboki, et al., "The P2Y(1) receptor antagonist adenosine-2',5'-diphosphate non-selectively antagonizes the platelet P2X(1) ion channel.", *Thromb. Haemostasis*, 86 (5), 1338 (2001)

Leonidas, et al., "Mapping the ribonucleolytic active site of eosinophil-derived neurotoxin (EDN) - High resolution crystal structures of EDN complexes with adenylic nucleotide inhibitors.", *J. Biol. Chem.*, 276 (18), 15009 (2001)

Gu, et al., "Crystal structures of the complexes of trichosanthin with four substrate analogs and catalytic mechanism of RNA N-glycosidase.", *Proteins*, 39 (1), 37 (2000)

1 unit = 1 µl at 10 mM solution

2'OMe-Adenosine-3',5'-bisphosphate

2'-O-Methyl-adenosine-3',5'-bisphosphate, Sodium salt

MW : 439.21 (Anion)

Description	Cat.#	Qty
2'OMe-Adenosine-3',5'-bisphosphate (2'OMe-pAp)	FQ6430	100 u

Guanosine-3',5'-bisphosphate

Guanosine-3',5'-bisphosphate, Triethylammonium salt

MW : 441.18 (Anion)

Description	Cat.#	Qty
Guanosine-3',5'-bisphosphate (pGp)	FQ2371	100 u
	FQ2370	500 u

Kvint, et al., "Emergency derepression: stringency allows RNA polymerase to override negative control by an active repressor.", *Mol. Microbiol.*, 35 (2), 435 (2000)

Acharya, et al., "The transmission of the electronic character of guanine-9-yl drives the sugar-phosphate backbone torsions in guanosine 3',5'-bisphosphate.", *Angew. Chem. Int. Edit.*, 38 (24), 3645 (1999)

Guanosine-2',5'-bisphosphate

Guanosine-2',5'-bisphosphate, Triethylammonium salt

MW: 441.18 (Anion)

Description	Cat.#	Qty
Guanosine-2',5'-bisphosphate	FQ2361	100 u
	FQ2360	500 u

Nucleic acid labeling and modification

Labeled and modified nucleotides

Xanthosine and Inosine nucleotides

Xanthosine nucleotides are used as G-nucleotide analogs in signal transduction research, molecular and cell biology.

Inosine is an analog of purine, and inosine nucleotides are widely used for mechanistic studies on ATP (cAMP) or GTP (cGMP) binding proteins.

XDP

Xanthosine-5'-diphosphate, Triethylammonium salt

MW : 442.17 (Anion)

Description	Cat.#	Qty
XDP	FQ2401	100 u
	FQ2400	500 u

Legate, *et al.*, "Nucleotide-dependent binding of the GTPase domain of the signal recognition particle receptor betasubunit to the alpha-subunit.", *J. Biol. Chem.*, 275 (35), 27439 (2000). Yu, *et al.*, "Interaction of the xanthine nucleotide binding Go alpha mutant with G protein-coupled receptors.", *J. Biol. Chem.* 273 (46), 30183 (1998)

XTP

Xanthosine-5'-triphosphate, Triethylammonium salt

MW : 521.14 (Anion)

Description	Cat.#	Qty
XTP	FQ2111	150 u
	FQ2110	750 u

Fulga, *et al.*, "SR beta coordinates signal sequence release from SRP with ribosome binding to the translocon.", *EMBO J.*, 20 (9), 2338(2001). Legate, *et al.*, "Nucleotide-dependent binding of the GTPase domain of the signal recognition particle receptor betasubunit to the alpha-subunit.", *J. Biol. Chem.* 275 (35), 27439 (2000)

IMP

Inosine-5'-monophosphate, Triethylammonium salt

MW : 347.20 (Anion)

Description	Cat.#	Qty
IMP	FQ2041	150 u
	FQ2040	750 u

Bhattacharya, *et al.*, "Regiospecificity of nucleotide-amino acid mating vs. water dynamics: a key to protein-nucleic acid assemblies: structure of unidecahydrated inosine-5'-monophosphate and L-glutamic acid.", *J. Chem. Crystallogr.* 30 (10), 655 (2000)
Heroux, *et al.*, "Crystal structures of the *Toxoplasma gondii* hypoxanthine-guanine phosphoribosyltransferase-GMP and -IMP complexes : Comparison of purine binding interactions with the XMP complex.", *Biochemistry-US* 38 (44), 14485 (1999)

IDP

Inosine-5'-diphosphate, Triethylammonium salt

MW : 426.17 (Anion)

Description	Cat.#	Qty
IDP	FQ2381	100 u
	FQ2380	500 u

AlAli, *et al.*, "Effects of metal ions on the activity of cytosolic phosphoenolpyruvate carboxykinase from camel kidney.", *Arab. Gulf J. Sci. Res.* 14 (3), 535 (1996). Vial, *et al.*, "Purification, partial kinetic characterization and reactive sulfhydryl-groups of the phosphoenolpyruvate carboxykinase from *peromyltilus-purpuratus* adductor muscle.", *Comp. Biochem. Phys. B.* 112 (3):451 (1995)

ITP

Inosine-5'-triphosphate, Triethylammonium salt

MW : 505.14

Description	Cat.#	Qty
ITP	FQ2051	150 u
	FQ2050	750 u

Noji, *et al.*, "Purine but not pyrimidine nucleotides support rotation of F(1)-ATPase.", *J. Biol. Chem.*, 276 (27), 25480 (2001)
Bianchi, *et al.*, "Intramolecular equilibria in metal ion complexes of guanosine 5'-triphosphate (GTP(4-)) and inosine 5'-triphosphate (ITP(4-)) in aqueous solution.", *J. Inorg. Biochem.*, 86 (1), 148 (2001). Chakrabarti, *et al.*, "Nucleoside triphosphate specificity of tubulin.", *Biochemistry*, 39 (33), 10269 (2000)

cIMP

Inosine-3',5'-cyclic monophosphate, Triethylammonium salt

MW : 329.18 (Anion)

Description	Cat.#	Qty
cIMP	FQ2850	50 u
	FQ2851	250 u

Wang, *et al.*, "RNA polymerase-cNMP-ligated cAMP receptor protein (CRP) mutant interactions in the enhancement of transcription by CRP mutants.", *J. Biol. Chem.*, 275 (43), 33457 (2000). Shapiro, *et al.*, "Structural basis for ligand selectivity of heteromeric olfactory cyclic nucleotide-gated channels.", *Biophys. J.*, 78 (5), 2307 (2000). Sunderman, *et al.*, "Sequence of events underlying the allosteric transition of rod cyclic nucleotide-gated channels.", *J. Gen. Physiol.*, 113 (5), 621 (1999)

1 u = 1 µl at 10 mM

Related products :

Description	Cat.#	Qty
XppNHp (XMPPNP)	FQ2990	50 u
XTPγS,	FQ3640	5 u
Mant-XDP	FP-FQ2941	10 u
Mant-XTP	FP-FQ2951	10 u
Mant-XppNHp (mant-XMPPNP)	FP-FQ5181	5 u
NPE-caged-XDP	FP-FQ3031	10 u
NPE-caged-XTP	FP-FQ3041	10 u
NPE-caged-XppNHp (NPE-caged-XMPPNP)	FP-FQ5201	5 u
lppNHp (lMPPNP)	FQ2890	50 u
Mant-EDA-ATP	FP-FQ5811	20 u
Mant-ITPgS	FP-FQ2931	10 u

D.148

interchim

Halogen-containing nucleotides

Halogen-containing nucleotides are used as rational phasing tools for protein crystallography, as specific substrates for studies on nucleotide binding proteins or as starting material for several nucleotide modifications.

5I-dUTP

5-Iodo-2'-deoxy-uridine-5'-triphosphate, Sodium salt
MW : 591.01 (Anion)

Description	Cat.#	Qty
5I-dUTP	FQ2600	50 u
	FQ2601	250 u

5I-dCTP

5-Iodo-2'-deoxy-cytidine-5'-triphosphate, Sodium salt
MW : 590.03 (Anion)

Description	Cat.#	Qty
5I-dCTP	FQ2590	50 u
	FQ2591	250 u

Bisubstrate inhibitors

Bisubstrate inhibitors are used for mechanistic studies on nucleotide kinases.

AP₃A

P¹-(5'-Adenosyl) P³-(5'-adenosyl) triphosphate, Sodium salt
MW : 753.38 (Anion)

Description	Cat.#	Qty
AP ₃ A	FQ2730	50 u
	FQ2731	250 u

Guranowski, *et al.*, "Selective degradation of 2'-adenylated diadenosine tri- and tetraphosphates, Ap(3)A and Ap(4)A, by two specific human dinucleoside polyphosphate hydrolases.", *Arch. Biochem. Biophys.*, 373, 218 (2000)

Luo, *et al.*, "Identification and characterization of diadenosine 5',5''-P₁,P₂-diphosphate and diadenosine 5',5''-P₁,P₃-triphosphate in human myocardial tissue.", *FASEB J.*, 13, 695 (1999)

Luthje, *et al.*, "Catabolism of Ap4A and Ap3A in whole blood. The dinucleotides are long-lived signal molecules in the blood ending up as intracellular ATP in the erythrocytes.", *Eur. J. Biochem.*, 173, 241 (1988)

AP₄A

P¹-(5'-Adenosyl) P⁴-(5'-adenosyl) tetraphosphate, Sodium salt
MW : 832.36 (Anion)

Description	Cat.#	Qty
AP ₄ A	FQ2301	100 u
	FQ2300	500 u

Vartanian, *et al.*, "Ap4A induces apoptosis in human cultured cells.", *FEBS Lett.*, 456, 175 (1999)

Campbell, *et al.*, "Characterization of P₁,P₄-diadenosine 5'-tetraphosphate binding on bovine aortic endothelial cells.", *Arch. Biochem. Biophys.*, 364, 280 (1999)

Guedon, *et al.*, "Effect of diadenosine tetraphosphate microinjection on heat shock protein synthesis in *Xenopus laevis* oocytes.", *EMBO J.*, 4, 3743 (1985)

AP₅A

P¹-(5'-Adenosyl) P⁵-(5'-adenosyl) pentaphosphate, Sodium salt
MW : 911.83 (Anion)

Description	Cat.#	Qty
AP ₅ A	FQ2760	50 u
	FQ2761	250 u

Wildman, *et al.*, "Selectivity of diadenosine polyphosphates for rat P2X receptor subunits.", *Eur. J. Pharmacol.*, 367, 119 (1999)

Pintor, *et al.*, "Presence of dinucleotide and ATP receptors in human cerebrocortical synaptic terminals.", *Eur. J. Pharmacol.* 366, 159 (1999)

Related products :

Description	Cat.#	Qty
5Br-dUDP	FQ4570	30 u
5Br-dUTP	FP-T8199A	25 µl

1 u = 1 µl at 10 mM

Nucleic acid labeling and modification

Labeled and modified nucleotides

AP₆A

P¹-(5'-Adenosyl) P⁶-(5'-adenosyl) hexaphosphate, Sodium salt

MW : 990.30 (Anion)

Description	Cat.#	Qty
AP ₆ A	FQ5600	20 u
	FQ5601	100 u

Morii, *et al.*, "Adenosine(5')hexaphospho(5')adenosine stimulation of a Ca(2+)-induced Ca(2+)-release channel from skeletal muscle sarcoplasmic reticulum.", *Eur. J. Biochem.*, 205, 979 (1992)

Luo, *et al.*, "Identification of diadenosine hexaphosphate in human erythrocytes.", *Hypertension*, 34, 872 (1999)

AP₅U

P¹-(5'-Adenosyl) P⁵-(5'-uridyl) pentaphosphate, Sodium salt

MW : 888.29 (Anion)

Description	Cat.#	Qty
AP ₅ U	FQ2790	50 u
	FQ2791	250 u

Cheng, *et al.*, "Homogeneous uridine kinase from Ehrlich ascites tumor: substrate specificity and inhibition by bisubstrate analogs.", *Mol Pharmacol.*, 30 (2), 159 (1986)

AP₄dT

P¹-(5'-adenosyl) P⁴-(5'-(2'-deoxy-thymidyl)) tetraphosphate, Triethylammonium salt

MW : 807.34 (Anion)

Description	Cat.#	Qty
AP ₄ dT	FQ2740	50 u
	FQ2741	250 u

AP₅dT

P¹-(5'-adenosyl) P⁵-(5'-(2'-deoxy-thymidyl)) pentaphosphate, Triethylammonium salt

MW : 886.32 (Anion)

Description	Cat.#	Qty
AP ₅ dT	FQ2770	50 u
	FQ2771	250 u

Lavie, *et al.*, "Structural basis for efficient phosphorylation of 3'-azidothymidine monophosphate by Escherichia coli thymidylate kinase.", *Proc. Natl. Acad. Sci., USA* 95 (24), 14045 (1998)

Lavie, *et al.*, "Crystal structure of yeast thymidylate kinase complexed with the bisubstrate inhibitor P-1-(5'-adenosyl) P-5-(5'-thymidyl) pentaphosphate (TP(5)A) at 2.0 angstrom resolution: Implications for catalysis and AZT activation.", *Biochemistry*, 37, 3677 (1998)

AP₄G

P¹-(5'-Adenosyl) P⁴-(5'-guanosyl) tetraphosphate, Triethylammonium salt

MW : 843.36 (Anion)

Description	Cat.#	Qty
AP ₄ G	FQ2750	50 u
	FQ2751	250 u

Ortiz, *et al.*, "specific synthesis of adenosine (5')tetraphospho(5')nucleoside and adenosine- (5')oligophospho(5')adenosine (n-greaterthan-4) catalyzed by firefly luciferase.", *Eur. J. Biochem.*, 212 (1), 263 (1993)

Palfi, *et al.*, "Alterations in the accumulation of adenylylated nucleotides in heavy-metal-ion-stressed and heat-stressed Synechococcus sp strain pcc-6301, a cyanobacterium, in light and dark.", *Biochem. J.* 276:487 (1991)

AP₅G

P¹-(5'-Adenosyl) P⁵-(5'-guanosyl) pentaphosphate, Triethylammonium salt

MW : 927.33 (Anion)

Description	Cat.#	Qty
AP ₅ G	FQ2780	50 u
	FQ2781	250 u

Prinz, *et al.*, "Binding of nucleotides to guanylate kinase, p21(ras), and nucleoside-diphosphate kinase studied by nano-electrospray mass spectrometry.", *J. Biol. Chem.*, 274 (50), 35337 (1999)

Ortiz, *et al.*, "specific synthesis of adenosine (5')tetraphospho(5')nucleoside and adenosine-(5')oligophospho(5')adenosine (n-greaterthan-4) catalyzed by firefly luciferase.", *Eur. J. Biochem.*, 212 (1), 263 (1993)

1 u = 1 µl at 10 mM

Antiviral nucleotides

AzTMP

3'-Azido-2',3'-dideoxy-thymidine-5'-monophosphate, Sodium salt

MW : 346.21 (Anion)

Description	Cat.#	Qty
AzTMP	FQ5620	20 u
	FQ5621	100 u

Cruchaga, *et al.*, "Inhibition of Phosphorolysis Catalyzed by HIV-1 Reverse Transcriptase Is Responsible for the Synergy Found in Combinations of 3'-Azido-3'-deoxythymidine with Nonnucleoside Inhibitors", *Biochemistry*, 44 (9), 3535 (2005)

Chenal-Francisque, *et al.*, "The highly similar TMP kinases of *Yersinia pestis* and *Escherichia coli* differ markedly in their AZTMP phosphorylating activity.", *Eur. J. Biochem.* 265, 112 (1999)

AzTTP

3'-Azido-2',3'-dideoxy-thymidine-5'-triphosphate, Sodium salt

MW : 504.16 (Anion)

Description	Cat.#	Qty
AzTTP	FQ2810	10 u
	FQ2811	50 u

Cruchaga, *et al.*, «Inhibition of Phosphorolysis Catalyzed by HIV-1 Reverse Transcriptase Is Responsible for the Synergy Found in Combinations of 3'-Azido-3'-deoxythymidine with Nonnucleoside Inhibitors», *Biochemistry*, 44 (9), 3535 (2005)

Akeb, *et al.*, «The production and evaluation of antibodies for enzyme immunoassay of AZTTP», *Nucleosides Nucleotides Nucleic Acids*, 20, 243 (2001) Faraj, *et al.*, «Effects of beta-L-3'-azido-3'-deoxythymidine 5'-triphosphate on host and viral DNA polymerases.», *Antiviral Res.* 47, 97 (2000)

d4TMP

2',3'-Didehydro-2',3'-dideoxy-thymidine-5'-monophosphate, Sodium salt

MW : 303.18 (Anion)

Description	Cat.#	Qty
d4TMP	FQ5640	20 u
	FQ5641	100 u

De Clercq, "Highlights in the development of new antiviral agents.", *Mini Rev. Med. Chem.*, 2, 163 (2002)

Mas, *et al.*, "Multidrug-resistant HIV-1 reverse transcriptase: involvement of ribonucleotide-dependent phosphorolysis in cross-resistance to nucleoside analogue inhibitors.", *J. Mol. Biol.*, 267, 181 (2002)

Selmi, *et al.*, "The valine-to-threonine 75 substitution in human immunodeficiency virus type 1 reverse transcriptase and its relation with stavudine resistance.", *J. Biol. Chem.*, 276, 13965 (2001)

d4TTP

2',3'-Didehydro-2',3'-dideoxy-thymidine-5'-triphosphate, Sodium salt

MW : 461.13 (Anion)

Description	Cat.#	Qty
d4TTP	FQ2861	10 u
	FQ2860	50 u

Ray, *et al.*, "Insights into the molecular mechanism of inhibition and drug resistance for HIV-1 RT with carbovir triphosphate.", *Biochemistry*, 41, 5150 (2002)

Vaccaro, *et al.*, "Mechanism of inhibition of the human immunodeficiency virus type 1 reverse transcriptase by d4TTP: an equivalent incorporation efficiency relative to the natural substrate dTTP.", *Antimicrob. Agents. Chemother.*, 44, 217 (2000)

3TCMP

2', 3'-Dideoxy-3'-thia-cytidine-5'-monophosphate, Sodium salt (L isomer)

MW : 308.22 (Anion)

Description	Cat.#	Qty
3TCMP	FQ2561	10 u
	FQ2560	50 u

Mas, *et al.*, "Multidrug-resistant HIV-1 reverse transcriptase: involvement of ribonucleotide-dependent phosphorolysis in cross-resistance to nucleoside analogue inhibitors.", *J. Mol. Biol.*, 323, 181 (2002)

3TCTP

2', 3'-Dideoxy-3'-thia-cytidine-5'-triphosphate, Sodium salt (L isomer)

MW : 466.16 (Anion)

Description	Cat.#	Qty
3TCTP	FQ3650	5 u
	FQ3651	25 u

Boyer, *et al.*, "YADD mutants of human immunodeficiency virus type 1 and Moloney murine leukemia virus reverse transcriptase are resistant to lamivudine triphosphate (3TCTP) in vitro.", *J. Virol.*, 75, 6321 (2001)

Kewn, *et al.*, "Development of enzymatic assays for quantification of intracellular lamivudine and carbovir triphosphate levels in peripheral blood mononuclear cells from human immunodeficiency virus-infected patients.", *Antimicrob. Agents Chemother.*, 46, 135 (2002)

1 u = 1 µl at 10 mM

Nucleic acid labeling and modification

Labeled and modified nucleotides

Ribavirin-5'-triphosphate

1-β-D-Ribofuranosyl-1,2,4-triazole-3-carboxamide-5'-triphosphate, Sodium salt
MW: 481.12 (Anion)

Description	Cat.#	Qty
Ribavirin-5'-triphosphate	FQ5940	20 u
	FQ5941	100 u

Maag, *et al.*, "Hepatitis C virus RNA-dependent RNA polymerase (NS5B) as a mediator of the antiviral activity of ribavirin.", *J. Biol. Chem.*, 276 (49), 46094 (2001)

Lanford, *et al.*, "Ribavirin induces error-prone replication of GB virus B in primary tamarin hepatocytes.", *J. Virol.*, 75 (17), 8074 (2001)

Crotty, *et al.*, "The broad-spectrum antiviral ribonucleoside ribavirin is an RNA virus mutagen.", *Nat. Med.*, 6 (12), 1375 (2000)

Other nucleotides analogs

N6-Methyl-ATP

N6-Methyl-adenosine-5'-triphosphate, Potassium salt
MW : 518.8 (Anion)

Description	Cat.#	Qty
N6-Methyl-ATP	FQ2101	150 u
	FQ2100	750 u

AP4

Adenosine-5'-tetrphosphate, Triethylammonium salt
MW : 583.13 (Anion)

Description	Cat.#	Qty
AP4	FQ2291	100 u
	FQ2290	500 u

ara-ATP

Adenine-arabinofuranoside-5'-triphosphate, Sodium salt
MW : 504.16 (Anion)

Description	Cat.#	Qty
ara-ATP	FQ2800	50 u
	FQ2801	250 u

Genini, *et al.*, "Nucleotide requirements for the in vitro activation of the apoptosis protein-activating factor-1-mediated caspase pathway.", *J. Biol. Chem.*, 275 (1), 29 (2000)

Ghoshal, *et al.*, "Ara-ATP impairs 3'-end processing of pre-messenger-RNAs by inhibiting both cleavage and polyadenylation.", *Nucleic Acids Res.* 19 (21):5871 (1991)

ε-ATP

1,N⁶-Etheno-adenosine-5'-triphosphate, Sodium salt
MW : 528.18 (Anion)

Description	Cat.#	Qty
ε-ATP	FQ5780	200 u
	FQ5781	100 u

Literature : Aguilar, *et al.*, "Ectoenzymatic breakdown of diadenosine polyphosphates by *Xenopus laevis* oocytes.", *Eur. J. Biochem.*, 268 (5), 1289 (2001)

Churchich, *et al.*, "A catalytic site of protein disulfide isomerase probed with adenosine-5'-triphosphate analogs.", *BBA-Protein Struct.M.*, 1479 (1-2), 293 (2000)

Gualix, *et al.*, "Studies of chromaffin granule functioning by flow cytometry : Transport of fluorescent epsilon-ATP and granular size increase induced by ATP.", *Receptor Channel*, 6 (6), 449 (1999)

8-Oxo-GTP

8-Oxo-guanosine-5'-triphosphate, Sodium salt
MW : 536.15 (Anion)

Description	Cat.#	Qty
8-Oxo-GTP	FQ2710	50 u
	FQ2711	250 u

Kuryavii, *et al.*, "Kinetics of inhibition of escherichia-coli RNA polymerase-catalyzed synthesis of the dinucleotide pppApU by 8-oxo-GTP and 8-BrGTP on the A1 promoter of Bacteriophage-17 delta-d111 DNA using a limited set of substrates.", *Mol. Biol. + 23* (3), 648 (1989)

1 u = 1 µl at 10 mM

8-Oxo-dGTP

8-Oxo-2'-deoxy-guanosine-5'-triphosphate, Sodium salt

MW : 519.15 (Anion)

Description	Cat.#	Qty
8-Oxo-dGTP	FQ4790	30 u
	FQ4791	150 u

Kai, et al., "An oxidized nucleotide affects DNA replication through activation of protein kinases in *Xenopus* egg lysates.", *Nucl. Acids Res.*, 30 (2), 569 (2002). Canitrot, et al., "Nucleotide excision repair DNA synthesis by excess DNA polymerase beta: a potential source of genetic instability in cancer cells.", *FASEB J.* 14 (12), 1765 (2000)
Nampalli, et al., "Efficient synthesis of 8-oxo-dGTP: A mutagenic nucleotide.", *Bioorg. Med. Chem. Lett.*, 10 (15), 1677 (2000)

6-Thio-GTP

6-Thio-guanosine-5'-triphosphate, Sodium salt

MW : 536.22 (Anion)

Description	Cat.#	Qty
6-Thio-GTP	FQ1991	150 u
	FQ1990	750 u

Tiede, et al., "CD28-dependent Rac1 activation is the molecular target of azathioprine in primary human CD4 + T lymphocytes.", *J. Clin. Invest.*, 111, 1133 (2003). Poland, et al., "Entrapment of 6-thiophosphoryl-IMP in the active site of crystalline adenylosuccinate synthetase from *Escherichia coli*.", *J. Biol. Chem.*, 272 (24), 15200 (1997)

5F-UTP

5-Fluoro-uridine-5'-triphosphate, Sodium salt

MW : 499.11 (Anion)

Description	Cat.#	Qty
5F-UTP	FQ2221	100 u
	FQ2220	500 u

Au, et al., "Reversed-phase ion-pair high-performance liquid-chromatographic assay of 5-fluorouracil, 5'-deoxy-5- fluorouridine, their nucleosides, monophosphate, diphosphate, and triphosphate nucleotides with a mixture of quaternary ammonium-ions.", *J. Chromatogr.*, 228, 245 (1982). Glazer, et al., "The effect of 5-fluorouridine 5'-triphosphate on RNA transcribed in isolated-nuclei in vitro.", *Mol.Pharmacol.* 17 (2), 279 (1980)

5-Aza-dCTP

5-Aza-2'-deoxy-cytidine-5'-triphosphate, Sodium salt

MW : 465.12 (Anion)

Description	Cat.#	Qty
5-Aza-dCTP	FQ5270	20 u
	FQ5271	100 u

6-Chlorpurine-riboside-5'-triphosphate

6-Chlorpurine-riboside-5'-triphosphate, Sodium salt

MW : 523.99 (Anion)

Description	Cat.#	Qty
6-Chlorpurine-riboside-5'-triphosphate	FQ2231	100 u
	FQ2230	500 u

Patzelt, et al., "Study of (Na+K)-ATPase with 6-chlorpurine triphosphate, 6-mercaptapurine triphosphate and dinitrophenyl-6-thiopurine triphosphate.", *H-S Z. Physiol. Chem.*, 355 (10), 1237 (1974)

6-Mercaptopurine-riboside-5'-triphosphate

6-Mercaptopurine-riboside-5'-triphosphate, Sodium salt

MW : 521.20 (Anion)

Description	Cat.#	Qty
6-Mercaptopurine-riboside-5'-triphosphate	FQ2241	100 u
	FQ2240	500 u

Patzelt, et al., "Study of (Na+K)-ATPase with 6-chlorpurine triphosphate, 6-mercaptapurine triphosphate and dinitrophenyl-6-thiopurine triphosphate.", *H-S Z. Physiol. Chem.*, 355 (10), 1237 (1974)

dPTP

6H,8H-3,4-Dihydro-pyrimido[4,5-c][1,2]oxazin-7-one-8-β-D-2'-deoxy-ribofuranosid-5'-triphosphate, Sodium salt

MW : 506.17 (Anion)

Description	Cat.#	Qty
dPTP	FQ4520	30 u
	FQ4521	150 u

Cummins, et al., "Synthesis and study of the fluorescein conjugate of the nucleotide dPTP.", *Nucleos. Nucleot.* 20 (4-7), 1049 (2001) ; Zaccolo, et al., "An approach to random mutagenesis of DNA using mixtures of Triphosphate derivatives of Nucleoside analogues.", *J. Mol. Biol.* 255 (4), 589 (1996)

1 u = 1 µl at 10 mM

Related products :

Description	Cat.#
ε-AppNHp (ε-AMPPNP)	FQ3660
Ribavirin-5'-triphosphate	FQ5940