

UptiTherm DNA Polymerase

Product Description

A superior polymerase to use as conventional polymerases but highly-thermostable, more accurate, and more...

UptiTherm Product format (buffer)	Qty	Cat. No. 1 U/μl	Cat. No. 5 U/μl
Standard buffer (including MgCl ₂)	1000 U	UPS53663	UPS53881
Mg free buffer + MgCl ₂ 50 mM	1000 U	UPS53703	UPS53921
Standard buffer (including MgCl ₂) + premixed dNTP (250 μl, 10 mM)	1000 U	UPS53733	UPS53763
Mg free buffer + MgCl ₂ 50 mM + premixed dNTP (250 μl, 10 mM)	1000 U	UPS53793	UPS53823

Note: UptiTherm is also available as 500 U and 250 U package on inquire, ans as a [gel form](#) (stable at room temperature!).

Store at - 20°C (1)

Technical and Scientific Information

Highly thermostable DNA polymerase. It is a recombinant, modified form of the enzyme from the thermophilic bacterium *Thermus thermophilus* expressed in *E. coli* (see Note 1).

UPTITHERM DNA polymerase is suitable for applications which require a highly thermostable and processive enzyme capable of synthesising DNA strands at elevated temperatures in amplification reactions or similar (e.g. primer extension), thus resolving the most complex secondary structures.

The enzyme is free of unspecific endo- or exonucleases activities, as well as nicking activities. It does not either exhibit significant reverse-transcriptase activity. Terminal transferase activity inherent to the enzyme renders A-tailed amplification products suitable to be further used in T/A cloning approaches.

The enzyme is supplied at a concentration of 5 U/μl and 1 U/μl in a storage buffer. 1 U/μl concentration allows accurate pipetting of small amounts of the DNA polymerase, so that it is not necessary to perform ulterior dilutions.

Unit Definition

One unit is defined as the amount of enzyme which incorporates 10 nanomoles of dNTPs into acid-insoluble DNA within 30 minutes at 72 °C.

Storage buffer

10 mM Tris-HCl (pH 8.0), 50 mM KCl, 1 mM EDTA, 0.1% Triton X-100, 50% glycerol (v/v). Other storage buffers available (contact our Technical Dpt.).

For any question,
contract your local distributor

Uptima,
powered by



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

uptima@interchim.com Hotline : 33(0)4 70 03 73 06

Reaction buffer

Recommended reaction buffer is: 75 mM Tris HCl (pH 9.0), 2 mM MgCl₂ (see Note 2), 50 mM KCl, 20 mM (NH₄)₂SO₄. This reaction buffer (the so-called Standard Buffer, Ordering Information at the end) is supplied at 10X concentration together with the enzyme.

Reaction buffer can be supplied MgCl₂ free (the so-called Free Buffer, see Ordering Information): Mg²⁺ ion, being the enzyme cofactor, plays a key role on polymerase activity, this is why its concentration must be optimized in certain amplification-based experiments. In this case, the MgCl₂ is supplied as a separate vial at 50 mM concentration. This solution must be completely thawed, vigorously vortexed and spun down in a bench-top centrifuge before use.

Storage Conditions

Store at -20°C in a constant temperature freezer (i.e. do not use frost-free freezers). Under these conditions the activity of the enzyme remains unaltered over 18 months of storage. The glycerol in the storage buffer prevents freezing at -20°C.

Reaction Conditions

After thawing the reaction buffer (and MgCl₂ solution, in case the “free buffer” choice is adopted), shake all vials (buffer, enzyme, 50 mM MgCl₂ solution) by gentle vortexing, later spin them down in a bench-top centrifuge, and eventually pipette desired volumes.

Keep all reagents on ice while they remain out of the -20°C storage freezer, otherwise enzyme activity will decrease over the time. Wear disposable gloves and make use of sterile, DNase- and RNase-free pipette tips and tubes in order to avoid contaminations and false negative results.

Recommended enzyme volumes to be added to the reaction mix

Final reaction volumes	Recommended enzyme volumes	Volume using 1U/μl items (or a 1U/μl master mix)
100 μl	up to 2.5U	Up to 2.5 μl
50 μl	1-1.25U	1-1.25 μl
25 μl	0.5-0.75U	0.5-0.75 μl

These volumes may vary depending on the reaction conditions and on the particular nucleotide composition of the product to be amplified. Further information may be obtained from our Technical Department.

It is recommended to increase the enzyme units (up to three times) in order to perform certain applications such as PRINS (Primed In Situ Synthesis) or when working on long DNA fragment amplifications (longer than 2 Kb from genomic DNA).

The dNTP final concentration recommended is 200 μM (50 μM each), but this figure may be decreased (e.g. when unspecific amplimers occur), enlarged (e.g. long amplifications) or even unbalanced in favour of any dNTP in particular (e.g. in vitro mutagenesis experiments) depending on the intended approach. This enzyme also accepts modified dNTPs (e.g. radioactively or fluorescein labelled) as substrate.

List of UptiTherm products

UPS53661	1000 units (1 U/μl) + Standard buffer (including 2 mM MgCl ₂)
UPS53662	500 units (1 U/μl) + Standard buffer (including 2 mM MgCl ₂)
UPS53663	250 units (1 U/μl) + Standard buffer (including 2 mM MgCl ₂)
UPS53703	1000 units (1 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 Mm
UPS53702	500 units (1 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM
UPS53701	250 units (1 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM
UPS53733	1000 units (1 U/μl) + Standard Buffer (with MgCl ₂) + premixed dNTP (250μl 10mM)
UPS53732	1000 units (1 U/μl) + Standard Buffer (with MgCl ₂) + premixed dNTP (10mM)
UPS53791	1000 units (1 U/μl) + Standard Buffer (with MgCl ₂) + premixed dNTP (10mM)
UPS53793	1000 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 Mm
UPS53792	500 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM
UPS53791	250 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM
UPS53881	1000 units (5 U/μl) + Standard buffer (including 2 mM MgCl ₂)
UPS53882	500 units (5 U/μl) + Standard buffer (including 2 mM MgCl ₂)
UPS53883	250 units (5 U/μl) + Standard buffer (including 2 mM MgCl ₂)
UPS53921	1000 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 Mm
UPS53922	500 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM
UPS53923	250 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM
UPS538763	1000 units (5 U/μl) + Standard Buffer (with MgCl ₂) + premixed dNTP (250μl 10mM)
UPS538762	1000 units (5 U/μl) + Standard Buffer (with MgCl ₂) + premixed dNTP (10mM)
UPS538761	1000 units (5 U/μl) + Standard Buffer (with MgCl ₂) + premixed dNTP (10mM)
UPS53823	1000 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 Mm
UPS53822	500 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM
UPS53821	250 units (5 U/μl) + Buffer MgCl ₂ free + vial MgCl ₂ 50 mM

Notes

Note 1: this enzyme is not recommended for certain experiments dealing with amplification of sequences homologous to those found in *E. coli* or very low-annealing temperature amplification approaches (e.g. RAPDs, Random Amplified Polymorphic DNAs). This enzyme is not recommended for sequencing.

Note 2: at difference with the vast majority of the thermostable DNA polymerases existing in the market, UPTHITHERM DNA polymerase shows optimal specificity at 2 mM MgCl₂ final concentration (rather than 1.5 mM) in reaction buffer.

Notice to buyers/users:

Some of the applications which may be performed with this product are covered by applicable patents in certain countries. The purchase of this product does not include or provide a license to perform patented applications. Users may be required to obtain a license depending on the country and/or application.

Rev.D0825E