




Advion Interchim
scientific



Bio-Works products

Research & Laboratory | Process Development | Bioprocess Production



Content

Introduction	4
Application areas	6
BabyBio Screening Kits	8
BabyBio™ NTA His-tag Screening Kit	
BabyBio IDA His-tag Screening Kit	
BabyBio IEX Screening Kit	
BabyBio Peptide Purification Kit	
Affinity chromatography	14
WorkBeads™ affimAb	
BabyBio affimAb	
WorkBeads Protein A	
BabyBio A	
Immobilized metal ion affinity chromatography	18
WorkBeads NiMAC and BabyBio NiMAC	
WorkBeads 40 NTA and WorkBeads 40 IDA precharged with Ni ²⁺ , Co ²⁺ , Zn ²⁺ , Cu ²⁺	
BabyBio NTA and BabyBio IDA precharged with Ni ²⁺ , Co ²⁺ , Zn ²⁺ , Cu ²⁺	
WorkBeads 40 NTA and WorkBeads 40 IDA uncharged	
BabyBio NTA and BabyBio IDA uncharged	
Ion exchange chromatography	30
Cation exchange	
WorkBeads 40S and WorkBeads 100S	
BabyBio S	
OptioBio 40S 10x100	
Anion exchange	
WorkBeads 40Q and WorkBeads 100Q	
BabyBio Q	
OptioBio™ 40Q 10x100	
WorkBeads 40 DEAE and BabyBio DEAE	
Multimodal ion exchange chromatography	38
WorkBeads 40 TREN	
BabyBio TREN	
Size exclusion chromatography	42
WorkBeads 40/100 SEC, WorkBeads 40/1000 SEC and	
WorkBeads 40/10 000 SEC	
WorkBeads Macro SEC	
WorkBeads 200 SEC	
Desalting/Buffer exchange	46
WorkBeads Dsalt and BabyBio Dsalt	
Pre-activated resins	50
WorkBeads 40/1000 ACT and WorkBeads 40/10 000 ACT	
BabyBio ACT	

WorkBeads, BabyBio and OptioBio are trademarks of Bio-Works Technologies AB. All other trademarks are the property of their respective owners. Products and technologies from Bio-Works Technologies AB are proprietary, patented or patent pending.





Introduction to Bio-Works

Uppsala, Sweden

Bio-Works is located at Uppsala Business Park site in Sweden. The company designs, develops, manufactures and supplies innovative leading edge products for purification and separation of peptides, proteins, oligonucleotides, viruses and other biomolecules for uses in research, process development and manufacturing in Life Science & Biopharma.

Bio-Works agarose based chromatography resins are manufactured using a proprietary method that results in porous beads with a tight size distribution and very high mechanical stability. Agarose based matrices have been successfully used for decades in biotechnology, from research to manufacturing scale purifications, due to its exceptional compatibility with biomolecules including proteins, peptides, and nucleic acids. Our company's WorkBeads resins are designed for separations that require optimal capacity, purity, productivity and reproducible scale-up results.

Experience and Quality

Bio-Works is highly experienced in the development and manufacturing of separation resins, and very knowledgeable about separation applications. The experience in our team covers the successful founding and growth of start-up companies, global and regional management within multinational companies, and the successful development, launch, and marketing of leading edge technologies.

The company is certified by Intertek and follows ISO 9001:2015.

Bio-Works supplies product information, quality documents, technical support, certificates, statements, vendor audits and regulatory support information.


Bio-Works believes in sustainability and care about the environment.




Technologies

Agar used as the starting material for WorkBeads resins is an inert, versatile and readily available natural material isolated from seaweed. It is the leading material used for purification matrices in protein science and biopharmaceutical processing. It will not denature or in any other way harm the delicate biotechnology products that are purified. Specialist knowledge is required to produce beads with suitable size, rigidity and porosity and then further to derivatize and make surface modifications for optimal final products.

The rigidity of the agarose based beads is important to avoid compression under high flow rates. Bio-Works patented cross-linking technology leverages high bead rigidity which allows very high flow rates. Large volumes can be processed fast and economically which is a key factor in manufacturing processes.





Products

Bio-Works advanced agarose based products are designed for purification of monoclonal and polyclonal antibodies, recombinant and native proteins, peptides, oligonucleotides, viruses, vaccines, enzymes, dairy proteins and for optimized purification of His-tagged proteins.

Bio-Works optimized WorkBeads resins are produced in several different bead sizes and porosities for both preparative research and bioprocess manufacturing scales. This allows seamless scalability and reproducible results. The bulk resins are available in pack sizes from 1.5 ml to 10 L and larger volumes on request.

The ready-to-use prepacked BabyBio columns (1 ml and 5 ml) and prepacked OptioBio 10x100 (7.9 ml) glass columns are designed for rapid, convenient and reproducible selectivity screenings and small scale purifications.

Several products are available for coupling of specific custom designed resins, for polishing of the target product in the final step, as well as, for very fast conditioning of the target product to prevent degradation.


Long term commitment

Bio-Works experience in agarose chemistry and long-term commitment ensures secured supply of products and continuous development of new chemistries, matrices and formats for future launches of high quality products for research, process development and manufacturing.

Our production and R&D departments are located in the same facility, this enables us to offer high flexibility and great technical service. In other words, we have the capacity and knowledge to develop and manufacture a large range of products optimized for many different application areas.

Bio-Works production meets your needs today and in the future. Our ambition is to make purification simple.





Application areas

Target molecule – how to start?



Antibodies

WorkBeads affimAb

- Top performance dynamic binding capacity also at short residence time
- Outstanding alkaline stability with 0.5 M NaOH
- WorkBeads 40 TREN for removing host cell proteins and extending the lifetime of the protein A resin

His-tagged proteins

Wide range of IMAC resins

- WorkBeads NiMAC with extra strongly bound Ni^{2+} resulting in extremely low nickel ion leakage, withstands 20 mM DTT and 20 mM EDTA
- Precharged and uncharged IMAC resins
- High selectivity, purity and yield
- NTA and IDA chelating ligands




Peptides

WorkBeads 40S and WorkBeads 40Q

- Higher dynamic binding capacity than other supplier's products
- Eliminates impurities prior to a polishing step, reducing bioburden on for example high-resolution silica resins
- BabyBio Peptide Purification Kit for fast and easy optimization





Oligonucleotides

WorkBeads 40Q

- Higher dynamic binding capacity compared to other supplier's products
- High purity when using standard low-pressure columns
- Applicable in large-scale industrial purification

Viruses and vaccines

WorkBeads 40/10 000 SEC


- WorkBeads 40S, WorkBeads 40Q and WorkBeads 40 TREN
- Purification of virus-like particles and virus components



Proteins & enzymes

WorkBeads 40S, WorkBeads 40Q and WorkBeads 40 TREN

- From capture to polishing
- Lower pressure with efficient purification
- High dynamic binding capacities at short residence times
- BabyBio IEX Screening Kit incl. four different 1 ml columns for fast and easy screening for running conditions and selectivity



BabyBio Screening Kits

BabyBio ready-to-use 1 ml and 5 ml columns are prepacked with most of Bio-Works' different resins. The column is made from biocompatible polypropylene which does not significantly interact with biomolecules. The top and bottom filters are made from polyethylene.

BabyBio Screening Kits are products for fast and easy screening of optimal selectivity and to find optimal running conditions in a reliable and convenient way.

Depending on purification scale when the optimal resin has been found for the specific target it is easy to continue using individual BabyBio columns or use bulk resin and pack a larger column.



BabyBio NTA His-tag Screening kits

Includes one 1 ml or one 5 ml of each: BabyBio Ni-NTA, BabyBio Co-NTA, BabyBio Cu-NTA and BabyBio Zn-NTA.

- Easy screening for optimized purity of His-tagged proteins
- Prepacked with WorkBeads 40 Ni-NTA, WorkBeads 40 Co-NTA, WorkBeads 40 Cu-NTA and WorkBeads 40 Zn-NTA

BabyBio IDA His-tag Screening kits

Includes one 1 ml or one 5 ml of each: BabyBio Ni-IDA, BabyBio Co-IDA, BabyBio Cu-IDA and BabyBio Zn-IDA.

- Easy screening for optimized purity of His-tagged proteins
- Prepacked with WorkBeads 40 Ni-IDA, WorkBeads 40 Co-IDA, WorkBeads 40 Cu-IDA and WorkBeads 40 Zn-IDA



BabyBio IEX Screening kit

Includes one 1 ml of each: BabyBio S, BabyBio Q, BabyBio DEAE and BabyBio TREN.


- Easy screening for optimal running conditions and selectivity
- Prepacked with WorkBeads 40S, WorkBeads 40Q, WorkBeads 40 DEAE and WorkBeads 40 TREN
- High binding capacity even at high flow rates
- Small-scale IEX purifications

BabyBio Peptide Purification Kit

Includes one 1 ml of each: BabyBio S and BabyBio Q

- Fast screening for optimal peptide purifications
- Prepacked with WorkBeads 40S and WorkBeads 40Q
- Easy screening for optimal running conditions
- High binding capacity even at high flow rates






Applications


Using BabyBio NTA His-tag Screening kit and BabyBio IDA His-tag Screening to find the optimal purification method

Sample: 2 ml His₆-GFP in binding buffer
Columns: BabyBio Ni-NTA, BabyBio Co-NTA, BabyBio Zn-NTA, BabyBio Cu-NTA
 BabyBio Ni-IDA, BabyBio Co-IDA, BabyBio Zn-IDA, BabyBio Cu-IDA
Column volume: 1 ml
Binding buffer: 50 mM sodium phosphate, 300 mM NaCl, 10 mM imidazole, pH 8.0
Elution buffer: 50 mM sodium phosphate, 300 mM NaCl, 300 mM imidazole, pH 8.0
Gradient: 0 to 100% elution buffer in 15 column volumes (CV)
Flow rate: 0.5 ml/min (75 cm/h)


BabyBio Ni-NTA 1 ml




BabyBio Ni-IDA 1 ml




BabyBio Co-NTA 1 ml




BabyBio Co-IDA 1 ml




BabyBio Cu-NTA 1 ml




BabyBio Cu-IDA 1 ml



BabyBio Zn-NTA 1 ml



BabyBio Zn-IDA 1 ml



Chromatograms showing comparisons of purifications of clarified His₆-GFP on BabyBio NTA 1 ml and BabyBio IDA 1 ml charged with Ni²⁺, Co²⁺, Cu²⁺ and Zn²⁺ ions. The blue and red lines correspond to the absorbance signal at 280 nm and 490 nm (specific for GFP), respectively, and the green line to the percentage of elution buffer.



Technical specifications

BabyBio NTA His-tag Screening kits and BabyBio IDA His-tag Screening kits

	BabyBio NTA His-tag Screening kit	BabyBio IDA His-tag Screening kit
Resin	WorkBeads 40 Ni-NTA, WorkBeads 40 Co-NTA WorkBeads 40 Cu-NTA, WorkBeads 40 Zn-NTA	WorkBeads 40 Ni-IDA, WorkBeads 40 Co-IDA WorkBeads 40 Cu-IDA, WorkBeads 40 Zn-IDA
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Ligand	Nitrilotriacetic acid (NTA)	Iminodiacetic acid (IDA)
Metal ion	Ni^{2+} , Co^{2+} , Cu^{2+} or Zn^{2+}	Ni^{2+} , Co^{2+} , Cu^{2+} or Zn^{2+}
Static binding capacity ²	70 mg His-tagged protein/ml resin	N/A
Dynamic binding capacity ² (DBC)	50 mg His-tagged protein/ml resin	N/A
Column volume	1 and 5 ml	1 and 5 ml
Column dimension	7 × 28 mm (1 ml), 13 × 38 mm (5 ml)	7 × 28 mm (1 ml), 13 × 38 mm (5 ml)
Recommended flow rates ³		
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ⁴		
BabyBio 1 ml	5 ml/min (780 cm/h)	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purifications, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 20% ethanol. Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped resin: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 100 mM sodium citrate-HCl (pH 3).	
pH stability	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² The binding capacity is determined using a BabyBio Ni-NTA 1 ml. The binding capacity is dependent on the size of the target protein, and on the competition with impurities.

³ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁴ Aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate at 4°C), or by additives (e.g., use half of the maximum flow rate for 20% ethanol).

BabyBio Peptide Purification kit

	BabyBio S	BabyBio Q
Target substance	Proteins, peptides	Protein, peptides, oligonucleotides, viruses
Resin	WorkBeads 40S	WorkBeads 40Q
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Ligand	Sulfonate ($-\text{SO}_3^-$)	Quaternary amine ($-\text{N}^+(\text{CH}_3)_3$)
Ion capacity	180 to 250 $\mu\text{mol H}^+$ /ml resin	180 to 250 $\mu\text{mol Cl}^-$ /ml resin
Dynamic binding capacity	130 mg BSA/ml resin ²	50 mg BSA/ml resin ³
Column volume	1 ml and 5 ml	1 ml and 5 ml
Column dimension	7 × 28 mm (1 ml), 13 × 38 mm (5 ml)	7 × 28 mm (1 ml), 13 × 38 mm (5 ml)
Recommended flow rate ⁴		
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ⁵		
BabyBio 1 ml	5 ml/min (780 cm/h)	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification. 1 M NaOH, 30% isopropanol and 70% ethanol. Should not be stored at low pH for prolonged time.	
pH stability	2 to 13	2 to 13
Storage	2 to 25°C in 20% ethanol with 0.2 M sodium acetate	2 to 25°C in 20% ethanol


¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 20 mM sodium citrate, 60 mM NaCl, pH 4.0.

³ Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 50 mM Tris-HCl, pH 8.0.

⁴ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁵ Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate for 20% ethanol).



BabyBio IEX Screening kit

	BabyBio S	BabyBio Q	BabyBio DEAE	BabyBio TREN
Target substance	Proteins, peptides	Protein, peptides, oligonucleotides, viruses	Protein, peptides, oligonucleotides	Proteins, peptides, oligonucleotides, viruses, chromatin fragments
Resin	WorkBeads 40S	WorkBeads 40Q	WorkBeads 40 DEAE	WorkBeads 40 TREN
Matrix	Rigid, highly cross-linked agarose			
Average particle size ¹ (D_{V50})	45 μ m	45 μ m	45 μ m	45 μ m
Ligand	Sulfonate (-SO ₃ ⁻)	Quaternary amine (-N ⁺ (CH ₃) ₃)	Diethylaminoethyl (-CH ₂ CH ₂ N ⁺ H(CH ₂ CH ₃) ₂)	Tris(2-aminoethyl) amine (TAEA)
Dynamic binding capacity	130 mg BSA/ml resin ²	50 mg BSA/ml resin ³	40 mg BSA/ml resin ³	50 mg BSA/ml resin ⁴
Column volume	1 ml and 5 ml	1 ml and 5 ml	1 ml and 5 ml	1 ml and 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rate ⁵				
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)	0.25 to 1 ml/min (37 to 150 cm/h)	0.25 to 1 ml/min (37 to 150 cm/h)	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)	1.25 to 5 ml/min (56 to 225 cm/h)	1.25 to 5 ml/min (56 to 225 cm/h)	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ⁶				
BabyBio 1 ml	5 ml/min (780 cm/h)	5 ml/min (780 cm/h)	5 ml/min (780 cm/h)	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)	20 ml/min (900 cm/h)	20 ml/min (900 cm/h)	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification. 1 M NaOH, 30% isopropanol and 70% ethanol. Should not be stored at low pH for prolonged time.			
pH stability	2 to 13	2 to 13	3 to 9 (rec. pH) 3 to 13	2 to 13
Storage	2 to 25°C in 20% ethanol with 0.2 M sodium acetate	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 20 mM sodium citrate, 60 mM NaCl, pH 4.0.

³ Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 50 mM Tris-HCl, pH 8.0.

⁴ Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 50 mM Tris-HCl, 50 mM NaCl, pH 8.0.

⁵ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁶ Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate for 20% ethanol).



Ordering information

Product name	Pack size	Article number
BabyBio NTA His-tag Screening kit 1 ml ¹	1 ml × 4	45 700 101
BabyBio NTA His-tag Screening kit 5 ml ¹	5 ml × 4	45 700 102
BabyBio IDA His-tag Screening kit 1 ml ¹	1 ml × 4	45 700 001
BabyBio IDA His-tag Screening kit 5 ml ¹	5 ml × 4	45 700 002
BabyBio IEX Screening Kit ²	1 ml × 4	45 900 001
BabyBio Peptide Purification Kit ³	1 ml × 2	45 300 102

¹ Includes one column each charged with Ni²⁺, Co²⁺, Cu²⁺ or Zn²⁺.

² Includes one 1 ml column of each: BabyBio S, BabyBio Q, BabyBio DEAE, BabyBio TREN.

³ Bundle of: BabyBio S 1 ml × 1 and BabyBio Q 1 ml × 1.

More information

Data Sheet, DS 40 650 010

WorkBeads 40 Ni-NTA, WorkBeads 40 Co-NTA,
WorkBeads 40 Zn-NTA, WorkBeads 40 Cu-NTA

WorkBeads 40 Ni-IDA, WorkBeads 40 Co-IDA,
WorkBeads 40 Zn-IDA, WorkBeads 40 Cu-IDA

BabyBio Ni-NTA, BabyBio Co-NTA,
BabyBio Zn-NTA, BabyBio Cu-NTA

BabyBio Ni-IDA, BabyBio Co-IDA,
BabyBio Zn-IDA, BabyBio Cu-IDA

BabyBio NTA His-tag Screening kit

BabyBio IDA His-tag Screening kit

www.bio-works.com/product/imac-resin

Data Sheet, DS 40 100 010


WorkBeads 40S, WorkBeads 40Q, WorkBeads 40 DEAE

BabyBio S, BabyBio Q, BabyBio DEAE, BabyBio TREN,
BabyBio IEX Screening Kit, BabyBio Peptide Purification Kit

www.bio-works.com/product/iex-resin








Affinity chromatography

Affinity chromatography (AC) separates proteins on the basis of a reversible interaction between a protein (or group of proteins) and a specific ligand coupled to a chromatography matrix. The technique is ideal for the capture step in a purification protocol. The target protein is collected in a highly pure and concentrated form. The high selectivity of affinity chromatography enables many purifications to be achieved in only one simple step, for example, purification of antibodies.

Target molecules

Monoclonal and polyclonal antibodies, bound via the Fc-region



Schematic depicting affinity chromatography


WorkBeads affimAb

- Top performance dynamic binding capacity also at short residence time
- Outstanding alkaline stability with 0.5 M NaOH, extends the number of purification cycles
- Excellent purity, recovery and reproducibility
- Negligible protein A leakage
- Convenient prepacked 1 ml and 5 ml BabyBio columns



WorkBeads Protein A


- For routine purification of antibodies in the research lab
- High dynamic binding capacity with excellent recovery and purity
- Reliable, reproducible and efficient
- Convenient prepacked 1 ml and 5 ml BabyBio columns



Applications

Dynamic binding capacity vs residence time

Resins: WorkBeads affimAb
MabSelect SuRe™ (Cytiva)
Column volume: 3.4 ml (6.6 × 100 mm)
Sample: 1 mg/ml human polyclonal IgG in PBS, pH 7.4
Binding buffer: PBS, pH 7.4
Elution buffer: 0.1 M glycine-HCl, pH 2.7
Cleaning-in-place (CIP): 5 column volumes (CV) 0.5 M NaOH at 2.4 min residence time (RT)
Residence times: 2.4, 4.8, 6 and 10 min (250, 125, 100 and 60 cm/h)



Alkaline stability comparison


Resins: WorkBeads affimAb
MabSelect SuRe
Column volume: 3.4 ml (6.6 × 100 mm)

DBC (10% breakthrough) determined at start and after each 20th CIP cycle

Sample: 1 mg/ml human polyclonal IgG in PBS, pH 7.4
Flow rate: 1.4 ml/min (2.4 min RT)
Binding buffer: PBS, pH 7.4
Elution buffer: 0.1 M glycine-HCl, pH 2.7



Each CIP cycle:

- 5 CV PBS, pH 7.4 at 1.4 ml/min (2.4 min RT)
- 0.5 M NaOH, 15 min contact time at 1 ml/min
- 5 CV PBS, pH 7.4 at 1.4 ml/min
- 5 CV 0.1 M glycine-HCl, pH 2.7 at 1.4 ml/min
- 5 CV PBS, pH 7.4 at 1.4 ml/min



HCP and HCD impurities in eluted mAb

Columns: WorkBeads affimAb: 6.6 × 50 mm, 1.7 ml
MabSelect SuRe: 6.6 × 50 mm, 1.7 ml
Feed: 18 ml clarified cell supernatant from CHO cells (100 cm/h)
Binding/wash buffer: PBS, pH 7.4 (300 cm/h)
Elution buffer: 100 mM glycine-HCl, pH 2.7 (150 cm/h)
CIP: 0.5 M NaOH (100 cm/h), 10 minutes contact time in each cycle
HCP: CHO CHP ELISA kit (#F550, Cygnus Technologies)
HCD: Quant-iT™ PicoGreen™ dsDNA Assay Kit (#P7589, ThermoFisher)





Technical specifications

	WorkBeads affimAb	WorkBeads Protein A
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	50 μm	45 μm
Ligand	Alkali stable recombinant protein A expressed in <i>E. coli</i> using animal-free medium	Recombinant protein A expressed in <i>E. coli</i> using animal-free medium
Dynamic binding capacity ² (DBC)	> 40 mg human IgG/ml resin	> 40 mg human IgG/ml resin
Max. recommended flow rate ³	300 cm/h	300 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 0.5 M NaOH, (pH 12), 10 mM HCl (pH 2), 0.1 M sodium citrate-HCl (pH 3), 6 M guanidine-HCl, 20% ethanol Should not be stored at low pH for prolonged time.	Compatible with all standard aqueous buffers used for protein purification, 10 mM HCl (pH 2), 0.1 M sodium citrate-HCl (pH 3), 6 M guanidine-HCl, 20% ethanol, 10 mM NaOH (pH 12) Should not be stored at low pH for prolonged time.
pH stability	3 to 10	3 to 10
Cleaning-in-place (CIP) stability	Up to 0.5 M NaOH	10 mM NaOH
Storage	2 to 8°C in 20 % ethanol	2 to 8°C in 20 % ethanol

¹ The median particle size of the cumulative volume distribution.

² DBC was determined at 10% breakthrough ($Q_{B10\%}$) by frontal analysis with 1 mg/ml human polyclonal IgG in PBS, pH 7.4 at 1.4 ml/min (245 cm/h, 2.5 minutes residence time) in a column packed with WorkBeads affimAb resin, column bed 6.6 × 100 mm.

³ Recommended flow rate at 20°C using aqueous buffers.


	BabyBio affimAb	BabyBio A
Resin	WorkBeads affimAb	WorkBeads Protein A
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	50 μm	45 μm
Ligand	Alkali stable recombinant protein A expressed in <i>E. coli</i> using animal-free medium	Recombinant protein A expressed in <i>E. coli</i> using animal-free medium
Dynamic binding capacity ² (DBC)	> 40 mg human IgG/ml resin	> 40 mg human IgG/ml resin
Column volume	1 ml 5 ml	1 ml 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rate ³		
BabyBio 1 ml	0.2 to 1 ml/min (28 to 150 cm/h)	0.3 to 1 ml/min (47 to 150 cm/h)
BabyBio 5 ml	0.9 to 4 ml/min (38 to 180 cm/h)	1 to 4 ml/min (45 to 180 cm/h)
Maximum flow rate ⁴		
BabyBio 1 ml	4 ml/min (620 cm/h)	4 ml/min (620 cm/h)
BabyBio 5 ml	15 ml/min (670 cm/h)	15 ml/min (670 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification	Compatible with all standard aqueous buffers used for protein purification
pH stability	3 to 10	3 to 10
Cleaning-in-place (CIP) stability	Up to 0.5 M NaOH	10 mM NaOH
Storage	2 to 8°C in 20 % ethanol	2 to 8°C in 20 % ethanol

¹ The median particle size of the cumulative volume distribution.

² DBC was determined at 10% breakthrough ($Q_{B10\%}$) by frontal analysis with 1 mg/ml human polyclonal IgG in PBS, pH 7.4 at 1.4 ml/min (245 cm/h, 2.5 minutes residence time) in a column packed with WorkBeads affimAb resin, column bed 6.6 × 100 mm.

³ Recommended flow rates include the flow rates in all steps; cleaning, equilibration, applying sample, washing, elution etc.

⁴ Decrease the max flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the max flow rate when operating at 4°C), or by additives (e.g., use half of the max flow rate for 20% ethanol).



Ordering information

Product name	Pack size	Article number
WorkBeads affimAb	25 ml	40 800 001
	200 ml	40 800 002
	1 L	40 800 010
	5 L	40 800 050
	10 L	40 800 060
BabyBio affimAb 1 ml	1 ml × 1	45 800 101
	1 ml × 2	45 800 102
	1 ml × 5	45 800 103
	1 ml × 10	45 800 104
BabyBio affimAb 5 ml	5 ml × 1	45 800 105
	5 ml × 2	45 800 106
	5 ml × 5	45 800 107
	5 ml × 10	45 800 108
WorkBeads Protein A	10 ml	40 605 003
	100 ml	40 605 004
	1 L	40 605 005
BabyBio A 1 ml	1 ml × 1	45 605 101
	1 ml × 2	45 605 102
	1 ml × 5	45 605 103
	1 ml × 10	45 605 104
BabyBio A 5 ml	5 ml × 1	45 605 105
	5 ml × 2	45 605 106
	5 ml × 5	45 605 107
	5 ml × 10	45 605 108

More information

Data Sheet, DS 40 800 010

WorkBeads affimAb, BabyBio affimAb

Data Sheet, DS 40 605 010

WorkBeads Protein A, BabyBio A

www.bio-works.com/product/affinity-chromatography




Immobilized metal ion affinity chromatography

Immobilized metal ion affinity chromatography (IMAC) separates most proteins with exposed histidine, cysteine and tryptophan on their surface. IMAC is an excellent technique for optimization and purification of His-tagged proteins. The technique is ideal for capture directly from clarified cell lysate. The target protein is collected in a highly purified and concentrated form.

Several factors influence the final purity of a His-tagged protein after an IMAC purification, for example, position of the tag (C- or N-terminal), the length of the tag, immobilized metal ion (Ni^{2+} , Co^{2+} , Zn^{2+} , Cu^{2+}) and the ligand immobilized on the matrix (NTA or IDA). To make the optimization of His-tagged protein purifications as efficient as possible Bio-Works offers products with many combinations of metal ion and immobilized ligand, as well as His-tag NTA Screening kits and His-tag IDA Screening kits.


Target molecules

His-tagged proteins and other proteins with exposed histidine, cysteine and tryptophan on their surface.



Schematic depicting Immobilized metal ion affinity chromatography





Precharged IMAC resin resistant to DTT and EDTA

WorkBeads NiMAC BabyBio NiMAC


- Resin with extra strongly bound Ni²⁺ resulting in extremely low nickel ion leakage
- Highly resistant to reducing agents up to 20 mM DTT
- Highly resistant to chelating substances present in eukaryotic extracts or up to 20 mM EDTA
- High purity and reproducible results
- Prepacked in BabyBio 1 ml and 5 ml columns




Applications

Larger sample load incl. 20 mM DTT and 20 mM Na₂-EDTA

Column: BabyBio NiMAC 1 ml
 Sample: 50 ml His₆-GFP in binding buffer with 20 mM DTT and 20 mM Na₂-EDTA
 Binding buffer: 50 mM sodium phosphate, 300 mM NaCl, 10 mM imidazole, pH 8.0
 Elution buffer: 50 mM sodium phosphate, 300 mM NaCl, 300 mM imidazole, pH 8.0
 Elution: Step gradient, 100% elution buffer, 10 column volumes (CV)
 Flow rates: 0.5 ml/min (78 cm/h; elution); 1 ml/min (loading)



(A) Chromatogram with 50 ml load of His₆-GFP, (B) SDS-PAGE under reducing conditions of the feed, flow-through (FT) and eluted pool (1x: concentrated eluate, 5x: 1:5 diluted eluate) from 10 ml sample load and 50 ml sample load. (C) Comparison of target recovery and purity for the two different sample load purifications.



IMAC resins precharged with different metal ions

WorkBeads 40 Ni-NTA
WorkBeads 40 Co-NTA
WorkBeads 40 Zn-NTA
WorkBeads 40 Cu-NTA
WorkBeads 40 Ni-IDA
WorkBeads 40 Co-IDA
WorkBeads 40 Zn-IDA
WorkBeads 40 Cu-IDA

- Resins immobilized with either NTA (Nitrilotriacetic acid) or IDA (Iminodiacetic acid) and four different choices of metal ions Ni^{2+} , Co^{2+} , Zn^{2+} or Cu^{2+}
- Precharged with different metal ions for ease of use
- Low leakage of immobilized ligand and metal ions
- Resistant to harsh cleaning agents (NaOH). *Note!* The metal ions have to be stripped off before cleaning
- High binding capacity and flow rate



BabyBio Ni-NTA, BabyBio Co-NTA
BabyBio Zn-NTA, BabyBio Cu-NTA


BabyBio Ni-IDA, BabyBio Co-IDA
BabyBio Zn-IDA, BabyBio Cu-IDA

- Prepacked for fast and reproducible purifications
- 1 ml and 5 ml columns with precharged WorkBeads NTA or WorkBeads IDA resins
- Easy to use with a syringe or chromatography system

BabyBio NTA His-tag Screening kit **BabyBio IDA His-tag Screening kit**

- Easy screening for optimized purity of His-tagged proteins
- 1 ml \times 4 and 5 ml \times 4 prepacked columns with precharged WorkBeads NTA or WorkBeads IDA resins for fast and convenient screening
- Each kit includes one column each of WorkBeads resins precharged with Ni^{2+} , Co^{2+} , Zn^{2+} and Cu^{2+}
- Easy to use with a syringe or chromatography system





Applications


Purification of His-tagged protein on BabyBio Ni-NTA vs. HisTrap FF +/- upstream BabyBio TREN

Purification of complex feeds expressed in different host cells can result in an extensive bioburden on the capture column in the form of DNA and different protein impurities when the feed is directly loaded without major pre-treatments. These impurities also often bind non-specifically to the target molecules and/or resin, and thus may be co-eluted with the final product.

One example of such a complex feed is His-tagged proteins expressed in bacteria. Therefore, WorkBeads 40 TREN prepacked in a BabyBio TREN column was introduced as a pre-treatment step upstream the IMAC columns when purifying the *E. coli* translation initiation factor 3 (IF-3-His₆). Since this protein has a nucleic acid binding domain, host cell nucleic acids can potentially be a major co-eluting impurity in the eluates.


The BabyBio TREN column was operated in flow-through mode to capture impurities. This resin binds host cell nucleic acid (HCD), endotoxins, viruses and various host cell proteins (HCP), thereby reducing the foulant load on the subsequent IMAC column. The effect of the TREN column was evaluated by examining the removal of HCD and HCPs, see analysis below.

Schematic view of BabyBio TREN upstream BabyBio Ni-NTA



Analytical SEC of eluted IF-3-His₆ from only BabyBio Ni-NTA compared to BabyBio Ni-NTA in combination with BabyBio TREN




Sample volumes: 200 μ L
 Column: Superdex™ 200 10/300 GL (Cytiva)
 Buffer: 20 mM PBS, pH 7.4
 Flow rate: 30 cm/h (0.3 ml/min)



IF-3-His₆ purification on BabyBio Ni-NTA and HisTrap FF +/- upstream BabyBio TREN


The effects of WorkBeads 40 TREN as a pre-treatment step:

- Increased purity in flow-throughs and eluates (Figures A-D)
- > 99.9% HCD removal in feed and 96-98% HCD removal in eluates
- No significant loss of target protein recovery (Figure B)
- ELISA showed that BabyBio TREN removes 49-62% more HCPs in eluates
- BabyBio Ni-NTA and HisTrap FF show analogous results (Figures B-D)



Conclusion

The final purity of the eluted His-tagged proteins, regarding HCD presence, was increased 97 – 98% when BabyBio TREN was placed upstream the IMAC column compared to a stand-alone step with only IMAC. Moreover, the purity was increased in terms of less co-eluted HCPs. In conclusion, the addition of an upstream pre-treatment step is highly advantageous to include early in a purification process to remove especially host-cell DNA that tend to interfere with downstream chromatographic processes.



Uncharged IMAC resins

WorkBeads 40 NTA

WorkBeads 40 IDA

- Resins immobilized with either NTA (Nitrilotriacetic acid) or IDA (Iminodiacetic acid) for immobilization of your choice of metal ion
- Low leakage of immobilized ligand and metal ions of choice
- Resistant to harsh cleaning agents (NaOH)
Note! The metal ions have to be stripped off before cleaning
- High binding capacity and flow rate



BabyBio NTA


BabyBio IDA

- Prepacked for fast and reproducible purifications
- 1 ml and 5 ml columns with WorkBeads NTA or WorkBeads IDA
- Fast and convenient to immobilize your choice of metal ions
- Easy to use with a syringe or chromatography system


Applications

Purification of clarified His₆-GFP on BabyBio Ni-NTA

Column: BabyBio Ni-NTA 1 ml
 Sample: 40 ml His₆-GFP in binding buffer
 Binding buffer: 50 mM sodium phosphate, 300 mM NaCl, 10 mM imidazole, pH 8.0
 Elution buffer: 50 mM sodium phosphate, 300 mM NaCl, 300 mM imidazole, pH 8.0
 Elution: 100% elution buffer in 5 CV
 Elution flow rate: 0.5 ml/min (75 cm/h)



(A) Chromatogram of the capture and elution of His₆-GFP. Absorbance at 280 nm (blue), absorbance at 490 nm (red) and percentage of elution buffer (green).
 (B) SDS-PAGE analysis of sample, flow-through, wash and eluted peak.



Technical specifications

WorkBeads NiMAC

Target substance	His-tagged proteins
Matrix	Highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm
Precharged ions	Nickel (II) ions, Ni^{2+}
Static binding capacity	> 80 mg/ml resin
Dynamic binding capacity ²	> 40 mg/ml resin
Metal ion capacity ³	> 60 $\mu\text{mol Cu}^{2+}$ /ml resin
Max flow rate (20 cm bed height and 5 bar)	600 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, and additives such as 20 mM $\text{Na}_2\text{-EDTA}$, 20 mM dithiothreitol (DTT), 20 mM TECP, 20 mM β -mercaptoethanol, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 500 mM imidazole, 30% isopropanol, 0.5 M NaOH
pH stability	3 to 9 (working range) 2 to 14 (cleaning-in-place)
Storage	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Binding capacity may vary depending on protein characteristics and on flow rate used. A lower flow rate usually increases the dynamic binding capacity.

³ Metal ion capacity is determined by frontal analysis at 50% breakthrough using copper solution.

BabyBio NiMAC

Target substance	His-tagged proteins
Matrix	WorkBeads NiMAC
Column volume	1 ml 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rate ¹	
BabyBio NiMAC 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio NiMAC 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)
Max flow rate ²	
BabyBio NiMAC 1 ml	5 ml/min (780 cm/h)
BabyBio NiMAC 5 ml	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, and additives such as 20 mM $\text{Na}_2\text{-EDTA}$, 20 mM dithiothreitol (DTT), 20 mM TECP, 20 mM β -mercaptoethanol, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 500 mM imidazole, 30% isopropanol, 0.5 M NaOH
pH stability	3 to 9 (working range) 2 to 14 (cleaning-in-place)
Storage	2 to 25°C in 20% ethanol

¹ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

² Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity.

Higher viscosities can be caused by low temperature or presence of additives. Use half of the maximum flow rate for 20% ethanol for example.



Technical specifications

	WorkBeads 40 Ni-NTA WorkBeads 40 Ni-IDA	WorkBeads 40 Co-NTA WorkBeads 40 Co-IDA	WorkBeads 40 Cu-NTA WorkBeads 40 Cu-IDA	WorkBeads 40 Zn-NTA WorkBeads 40 Zn-IDA
Matrix	Highly cross-linked agarose	Highly cross-linked agarose	Highly cross-linked agarose	Highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm	45 μm	45 μm
Chelating ligand	Nitrilotriacetic acid (NTA) or Iminodiacetic acid (IDA)	NTA or IDA	NTA or IDA	NTA or IDA
Metal ion	Nickel (II)	Cobalt (II)	Copper (II)	Zink (II)
Metal ion capacity for the chelating ligand ²	N/A	N/A	50 to 60 $\mu\text{mol Cu}^{2+}/\text{ml}$ (WorkBeads 40 Cu-IDA)	N/A
Dynamic binding capacity ³ (DBC)	> 60 mg His ₆ -GFP/ml resin	N/A	N/A	N/A
Maximum flow rate (20 cm bed height, 5 bar)	600 cm/h	600 cm/h	600 cm/h	600 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 20% ethanol. Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped resin: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 10 mM sodium citrate-HCl (pH 3).			
pH stability	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Metal ion capacity is determined by frontal analysis at 50% breakthrough using copper solution.

³ The binding capacity is determined using a BabyBio Ni-NTA 1 ml, equal value is expected for IDA resins. The binding capacity is dependent on the size of the target protein, and on the competition of impurities.


	BabyBio: Ni-NTA, Co-NTA, Cu-NTA, Zn-NTA	BabyBio: Ni-IDA, Co-IDA, Cu-IDA, Zn-IDA
Resin	WorkBeads 40 Ni-NTA, WorkBeads 40 Co-NTA WorkBeads 40 Cu-NTA, WorkBeads 40 Zn-NTA	WorkBeads 40 Ni-IDA, WorkBeads 40 Co-IDA WorkBeads 40 Cu-IDA, WorkBeads 40 Zn-IDA
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Ligand	Nitrilotriacetic acid (NTA)	Iminodiacetic acid (IDA)
Metal ion	Ni ²⁺ , Co ²⁺ , Cu ²⁺ or Zn ²⁺	Ni ²⁺ , Co ²⁺ , Cu ²⁺ or Zn ²⁺
Static binding capacity ²	70 mg His-tagged protein/ml resin	N/A
Dynamic binding capacity ² (DBC)	50 mg His-tagged protein/ml resin	N/A
Column volume	1 and 5 ml	1 and 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rates ³		
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ⁴		
BabyBio 1 ml	5 ml/min (780 cm/h)	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purifications, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 20% ethanol. Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped resin: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 100 mM sodium citrate-HCl (pH 3).	
pH stability	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² The binding capacity is determined using a BabyBio Ni-NTA 1 ml. The binding capacity is dependent on the size of the target protein, and on the competition with impurities.

³ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁴ Aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate at 4°C), or by additives (e.g., use half of the maximum flow rate for 20% ethanol).



Technical specifications

	WorkBeads 40 NTA	WorkBeads 40 IDA
Matrix	Highly cross-linked agarose	Highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Chelating ligand	Nitrilotriacetic acid (NTA)	Iminodiacetic acid (IDA)
Metal ion capacity ²	20 to 30 $\mu\text{mol Cu}^{2+}/\text{ml resin}$	50 to 60 $\mu\text{mol Cu}^{2+}/\text{ml resin}$
Maximum flow rate	600 cm/h (20 cm bed height, 5 bar)	600 cm/h (20 cm bed height, 5 bar)
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 20% ethanol. Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped resin: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 10 mM sodium citrate-HCl (pH 3).	
pH stability	7 to 9 (working range) 2 to 12 (cleaning, stripped resin)	7 to 9 (working range) 2 to 12 (cleaning, stripped resin)
Storage	2 to 25°C in ethanol	2 to 25°C in ethanol

¹ The median particle size of the cumulative volume distribution.

² Metal ion capacity is determined by frontal analysis at 50% breakthrough using copper solution.

	BabyBio NTA	BabyBio IDA
Resin	WorkBeads 40 NTA	WorkBeads 40 IDA
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Ligand	NTA	IDA
Column volume	1 ml 5 ml	1 ml 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rates ²		
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ³		
BabyBio 1 ml	5 ml/min (780 cm/h)	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purifications, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 20% ethanol. Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped resin: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 100 mM sodium citrate-HCl (pH 3).	
pH stability	7 to 9 (working range) 2 to 12 (cleaning, stripped resin)	7 to 9 (working range) 2 to 12 (cleaning, stripped resin)
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.


³ Aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate at 4°C), or by additives (e.g. use half of the maximum flow rate for 20% ethanol).



Ordering information

Product name	Pack size	Article number
WorkBeads NiMAC	25 ml	40 653 001
	150 ml	40 653 003
	1 L	40 653 010
WorkBeads 40 Ni-NTA	25 ml	40 651 001
	150 ml	40 651 003
	1 L	40 651 010
WorkBeads 40 Co-NTA	25 ml	40 651 401
	150 ml	40 651 403
	1 L	40 651 410
WorkBeads 40 Cu-NTA	25 ml	40 651 301
	150 ml	40 651 303
	1 L	40 651 310
WorkBeads 40 Zn-NTA	25 ml	40 651 501
	150 ml	40 651 503
	1 L	40 651 510
WorkBeads 40 Ni-IDA	25 ml	40 650 001
	150 ml	40 650 003
	1 L	40 650 010
WorkBeads 40 Co-IDA	25 ml	40 650 401
	150 ml	40 650 403
	1 L	40 650 410
WorkBeads 40 Cu-IDA	25 ml	40 650 301
	150 ml	40 650 303
	1 L	40 650 310
WorkBeads 40 Zn-IDA	25 ml	40 650 501
	150 ml	40 650 503
	1 L	40 650 510
BabyBio NiMAC 1 ml	1 ml × 1	45 655 311
	1 ml × 2	45 655 312
	1 ml × 5	45 655 313
	1 ml × 10	45 655 314
BabyBio NiMAC 5 ml	5 ml × 1	45 655 315
	5 ml × 2	45 655 316
	5 ml × 5	45 655 317
	5 ml × 10	45 655 318
BabyBio NTA His-tag Screening kit 1 ml ¹	1 ml × 4	45 700 101
BabyBio NTA His-tag Screening kit 5 ml ¹	5 ml × 4	45 700 102
BabyBio IDA His-tag Screening kit 1 ml ¹	1 ml × 4	45 700 001
BabyBio IDA His-tag Screening kit 5 ml ¹	5 ml × 4	45 700 002
BabyBio Ni-NTA 1 ml	1 ml × 1	45 655 101
	1 ml × 2	45 655 102
	1 ml × 5	45 655 103
	1 ml × 10	45 655 104
	1 ml × 100	45 655 110
BabyBio Ni-NTA 5 ml	5 ml × 1	45 655 105
	5 ml × 2	45 655 106
	5 ml × 5	45 655 107
	5 ml × 10	45 655 108
	5 ml × 100	45 655 109
BabyBio Co-NTA 1 ml	1 ml × 1	45 655 131
	1 ml × 2	45 655 132
	1 ml × 5	45 655 133
	1 ml × 10	45 655 134
BabyBio Co-NTA 5 ml	5 ml × 1	45 655 135
	5 ml × 2	45 655 136
	5 ml × 5	45 655 137
	5 ml × 10	45 655 138
BabyBio Cu-NTA 1 ml	1 ml × 1	45 655 121
	1 ml × 2	45 655 122
	1 ml × 5	45 655 123
	1 ml × 10	45 655 124
BabyBio Cu-NTA 5 ml	5 ml × 1	45 655 125
	5 ml × 2	45 655 126
	5 ml × 5	45 655 127
	5 ml × 10	45 655 128

¹ Includes one column each charged with Ni²⁺, Co²⁺, Cu²⁺ or Zn²⁺



Ordering information

Product name	Pack size	Article number
BabyBio Zn-NTA 1 ml	1 ml × 1	45 655 141
	1 ml × 2	45 655 142
	1 ml × 5	45 655 143
	1 ml × 10	45 655 144
BabyBio Zn-NTA 5 ml	5 ml × 1	45 655 145
	5 ml × 2	45 655 146
	5 ml × 5	45 655 147
	5 ml × 10	45 655 148
BabyBio Ni-IDA 1 ml	1 ml × 1	45 655 001
	1 ml × 2	45 655 002
	1 ml × 5	45 655 003
	1 ml × 10	45 655 004
BabyBio Ni-IDA 5 ml	5 ml × 1	45 655 005
	5 ml × 2	45 655 006
	5 ml × 5	45 655 007
	5 ml × 10	45 655 008
BabyBio Co-IDA 1 ml	1 ml × 1	45 655 031
	1 ml × 2	45 655 032
	1 ml × 5	45 655 033
	1 ml × 10	45 655 034
BabyBio Co-IDA 5 ml	5 ml × 1	45 655 035
	5 ml × 2	45 655 036
	5 ml × 5	45 655 037
	5 ml × 10	45 655 038
BabyBio Cu-IDA 1 ml	1 ml × 1	45 655 021
	1 ml × 2	45 655 022
	1 ml × 5	45 655 023
	1 ml × 10	45 655 024
BabyBio Cu-IDA 5 ml	5 ml × 1	45 655 025
	5 ml × 2	45 655 026
	5 ml × 5	45 655 027
	5 ml × 10	45 655 028
BabyBio Zn-IDA 1 ml	1 ml × 1	45 655 041
	1 ml × 2	45 655 042
	1 ml × 5	45 655 043
	1 ml × 10	45 655 044
BabyBio Zn-IDA 5 ml	5 ml × 1	45 655 045
	5 ml × 2	45 655 046
	5 ml × 5	45 655 047
	5 ml × 10	45 655 048

More information

Data Sheet, DS 40 653 010

WorkBeads NiMAC, BabyBio NiMAC

Data Sheet, DS 40 650 010

WorkBeads 40 Ni-NTA, WorkBeads 40 Co-NTA,
WorkBeads 40 Zn-NTA, WorkBeads 40 Cu-NTA

WorkBeads 40 Ni-IDA, WorkBeads 40 Co-IDA,
WorkBeads 40 Zn-IDA, WorkBeads 40 Cu-IDA

BabyBio Ni-NTA, BabyBio Co-NTA,
BabyBio Zn-NTA, BabyBio Cu-NTA
BabyBio Ni-IDA, BabyBio Co-IDA,
BabyBio Zn-IDA, BabyBio Cu-IDA

BabyBio NTA His-tag Screening kit
BabyBio IDA His-tag Screening kit

www.bio-works.com/product/imac-resin





Ordering information

Product name	Pack size	Article number
WorkBeads 40 NTA	25 ml	40 602 001
	150 ml	40 602 003
	1 L	40 602 010
WorkBeads 40 IDA	25 ml	40 601 001
	150 ml	40 601 003
	1 L	40 601 010
BabyBio NTA 1 ml	1 ml × 1	45 655 111
	1 ml × 2	45 655 112
	1 ml × 5	45 655 113
	1 ml × 10	45 655 114
BabyBio NTA 5 ml	5 ml × 1	45 655 115
	5 ml × 2	45 655 116
	5 ml × 5	45 655 117
	5 ml × 10	45 655 118
BabyBio IDA 1 ml	1 ml × 1	45 655 011
	1 ml × 2	45 655 012
	1 ml × 5	45 655 013
	1 ml × 10	45 655 014
BabyBio IDA 5 ml	5 ml × 1	45 655 015
	5 ml × 2	45 655 016
	5 ml × 5	45 655 017
	5 ml × 10	45 655 018

More information


Data Sheet, DS 40 600 010


WorkBeads 40 NTA, WorkBeads 40 IDA

BabyBio NTA, BabyBio IDA

www.bio-works.com/product/imac-resin








Ion exchange chromatography

Ion exchange chromatography (IEX) is a very useful technique that separates proteins on the basis of differences in their net surface charge in relation to pH of the surroundings. Every protein has its own charge/pH relationship.

WorkBeads ion exchangers are ideal for proteins, peptides and oligonucleotides. They show excellent results for larger peptides, particularly insulin. Two different bead sizes are available for optimal purity in all different steps during a purification process, the capture, the enhancement and the polishing steps.

Target molecules

In general most proteins, peptides and oligonucleotides. IEX is a universal purification technique suitable in all purification steps in a process.




Schematic depicting ion exchange chromatography




Cation exchangers

WorkBeads 40S WorkBeads 100S

- Designed for research and industrial scale purifications of proteins and peptides
- Two different bead sizes, 45 μm and 100 μm
- High chemical stability for easy cleaning-in-place
- High binding capacity during high flow rates



Structure of the ligand used in WorkBeads 40S



OptioBio 40S 10x100

- Prepacked glass column for reliable and reproducible results
- Optimal for high-performance small-scale purification
- 10 cm bed height for fast method optimization in bioprocess development



BabyBio S

- Prepacked for fast and reproducible purifications
- BabyBio 1 ml and 5 ml columns
- Easy to use with a syringe or chromatography system


Applications

Comparison of prepacked cation exchangers


Prepacked column: OptioBio 40S 10x100
 Resin: WorkBeads 40S
 Sample: 2.5 ml 1.5 mg/ml Concanavalin A, 1.5 mg/ml Ribonuclease A, 0.5 mg/min α -chymotrypsinogen A, 0.5 mg/ml Lysozyme
 Binding buffer: 50 mM MES, pH 6.0
 Elution buffer: 0-50% in 20 CV 50 mM MES, 1 M NaCl, pH 6.0
 Flow rate: 2 ml/min, 150 cm/h, 4 min RT

Prepacked column: HiScreen™ Capto™ SP ImpRes (Cytiva)
 Resin: Capto SP ImpRes
 Sample: 1.5 ml 1.5 mg/ml Concanavalin A, 1.5 mg/ml Ribonuclease A, 0.5 mg/min α -chymotrypsinogen A, 0.5 mg/ml Lysozyme
 Binding buffer: 50 mM MES, pH 6.0
 Elution buffer: 0-50% in 20 CV 50 mM MES, 1 M NaCl, pH 6.0
 Flow rate: 1.2 ml/min, 150 cm/h, 4 min RT

OptioBio 40S 10x100



HiScreen Capto SP ImpRes



Peptide purification, comparison of dynamic binding capacity and purity

Sample: 45 amino acid residue peptide
 Resins: WorkBeads 40S
 Capto SP ImpRes
 Column: 10 × 240 mm, 19 ml
 Flow: 2 ml/min (150 cm/h)
 Buffers: 15% acetonitrile in a proprietary buffer composition

Resin	DBC (mg/ml) at 2.0 min residence time	DBC (mg/ml) at 1.1 min residence time	Purity (%) ¹
WorkBeads 40S	150	140	91.8
Capto SP ImpRes	125	123	85.2

¹ Load of 30 g/L crude feed containing 55% target peptide.



Technical specifications

	WorkBeads 40S	WorkBeads 100S
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	90 to 110 μm
Ionic group (ligand)	Sulfonate ($-\text{SO}_3^-$)	Sulfonate ($-\text{SO}_3^-$)
Ionic capacity	180 to 250 $\mu\text{mol H}^+/\text{ml resin}$	180 to 250 $\mu\text{mol H}^+/\text{ml resin}$
Dynamic binding capacity ² (DBC)	130 mg BSA/ml resin	> 100 mg BSA/ml resin
Pressure flow characteristic	N/A	2 bar at 900 cm/h, 25 mm diameter \times 20 cm bed height
Maximum flow rate	600 cm/h (20 cm bed height, 5 bar)	N/A
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 1 M NaOH, 30% isopropanol or 70% ethanol. Should not be stored at < pH 3 for prolonged time.	
pH stability	2 to 13	2 to 13
Storage	2 to 25°C in 20% ethanol with 0.2 M sodium acetate	2 to 25°C in 20% ethanol with 0.2 M sodium acetate

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4-minutes residence time in 20 mM Na-citrate, pH 4.0.


OptioBio 40S 10x100

Resin	WorkBeads 40S
Matrix	Rigid, highly cross-linked agarose
Average particle size (D_{v50}) ¹	45 μm
Ionic group (ligand)	Sulfonate ($-\text{SO}_3^-$)
Ionic capacity	180 to 250 $\mu\text{mol H}^+/\text{ml resin}$
Dynamic binding capacity ² (DBC)	150 mg BSA/ml resin
Column volume	7.9 ml
Column dimension	10 \times 100 mm
Recommended flow rate	2 to 4 ml/min (150 to 300 cm/h)
Maximum flow rate ³	6 ml/min (450 cm/h)
Column hardware pressure limit	2.1 MPa, 21 bar, 305 psi
Chemical stability	Compatible with all standard buffers used for protein purification, 1 M NaOH, 30 % isopropanol or 70 % ethanol. Should not be stored at < pH 3 for prolonged time.
pH stability	2 to 13
Storage	2 to 25°C in 20% ethanol with 0.2 M sodium acetate

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined in 20 mM Na-citrate, pH 4.0, at a flow of 2 ml/min (150 cm/h; 4 minutes residence time).

³ Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature. Use half of the maximum flow rate for 20% ethanol.



BabyBio S	
Resin	WorkBeads 40S
Matrix	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm
Ionic group (ligand)	Sulfonate ($-\text{SO}_3^-$)
Ion capacity	180 to 250 $\mu\text{mol H}^+/\text{ml resin}$
Dynamic binding capacity ² (DBC)	130 mg BSA/ml resin
Column volume	1 ml and 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rate ³	
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)

BabyBio S	
Maximum flow rate ⁴	
BabyBio 1 ml	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification and 70% ethanol.
	Should not be stored at low pH for prolonged time.
pH stability	2 to 13
Storage	2 to 25°C in 20% ethanol with 0.2 M sodium acetate

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 20 mM Na-citrate, 60 mM NaCl, pH 4.0.

³ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁴ Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate for 20% ethanol).

Ordering information

Product name	Pack size	Article number
WorkBeads 40S	25 ml	40 200 001
	200 ml	40 200 002
	1 L	40 200 010
	5 L	40 200 050
	10 L	40 200 060
WorkBeads 100S	25 ml	10 200 001
	200 ml	10 200 002
	500 ml	10 200 005
	1 L	10 200 010
	5 L	10 200 050
	10 L	10 200 060
OptioBio 40S 10×100	7.9 ml × 1	55 420 011

Product name	Pack size	Article number
BabyBio IEX Screening Kit ¹	1 ml × 4	45 900 001
BabyBio Peptide Purification Kit ²	1 ml × 2	45 300 102
BabyBio S 1 ml	1 ml × 1	45 200 101
	1 ml × 2	45 200 102
	1 ml × 5	45 200 103
	1 ml × 10	45 200 104
BabyBio S 5 ml	5 ml × 1	45 200 105
	5 ml × 2	45 200 106
	5 ml × 5	45 200 107
	5 ml × 10	45 200 108

¹ Includes one 1 ml column of each: BabyBio S, BabyBio Q, BabyBio DEAE, BabyBio TREN.

² Bundle of: BabyBio S 1 ml × 1 and BabyBio Q 1 ml × 1.

More information

Data Sheet, DS 40 100 010

WorkBeads 40S, WorkBeads 40Q,
WorkBeads 40 DEAE

BabyBio S, BabyBio Q, BabyBio DEAE,
BabyBio TREN, BabyBio IEX Screening Kit,
BabyBio Peptide Purification Kit


www.bio-works.com/product/iex-resin
www.bio-works.com/product/optiobio-columns

Data Sheet, DS 10 200 010

WorkBeads 100S, WorkBeads 100Q

Data Sheet, DS 55 410 010

OptioBio 40S 10×100, OptioBio 40Q 10×100





Anion exchangers



Target molecules

In general most proteins and peptides. IEX is a universal purification technique suitable in all purification steps in a process. Anion exchangers are also excellent for purification of oligonucleotides.



WorkBeads 40Q WorkBeads 100Q WorkBeads 40 DEAE

- Available with two different ligands, strong anion exchanger (Q) and weak anion exchanger (DEAE)
- Designed for research and industrial scale purifications of proteins, peptides and oligonucleotides
- Two different bead sizes, 45 μm and 100 μm
- High chemical stability for easy cleaning-in-place
- High binding capacity during high flow rates



OptioBio 40Q 10x100


- Prepacked glass column for reliable and reproducible results
- Optimal for high-performance small-scale purification
- 10 cm bed height for fast method optimization in bioprocess development



BabyBio Q BabyBio DEAE

- Prepacked for fast and reproducible purifications
- 1 ml and 5 ml columns packed with WorkBeads 40Q and WorkBeads 40 DEAE
- Easy to use with a syringe or chromatography system






Applications

Comparison of prepacked anion exchangers


Prepacked column: OptioBio 40Q 10x100
 Resin: WorkBeads 40Q
 Sample: 10 ml 0.7 mg/ml apo-transferrin, 0.45 mg/ml α -lactalbumin, 1.4 mg/ml soybean trypsin inhibitor
 Binding buffer: 50 mM Tris-HCl, pH 7.4
 Elution buffer: 0-40% over 20 CV, 50 mM Tris-HCl, 1 M NaCl pH 7.4
 Flow rate: 2 ml/min, 150 cm/h, 4 min RT

Prepacked column: HiScreen Capto Q ImpRes (Cytiva)
 Resin: Capto Q ImpRes
 Sample: 6 ml 0.7 mg/ml apo-transferrin, 0.45 mg/ml α -lactalbumin, 1.4 mg/ml soybean trypsin inhibitor
 Binding buffer: 50 mM Tris-HCl, pH 7.4
 Elution buffer: 0-40% over 20 CV, 50 mM Tris-HCl, 1 M NaCl pH 7.4
 Flow rate: 1.2 ml/min, 150 cm/h, 4 min RT

OptioBio 40Q 10x100



HiScreen Capto Q ImpRes



Technical specifications

	WorkBeads 40Q	WorkBeads 100Q
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μ m	90 to 110 μ m
Ionic group (ligand)	Quarternary amine ($-N^+(CH_3)_3$)	Quarternary amine ($-N^+(CH_3)_3$)
Ionic capacity	180 to 250 μ mol Cl ⁻ /ml resin	140 to 200 μ mol Cl ⁻ /ml resin
Dynamic binding capacity ² (DBC)	47 mg BSA/ml resin	> 40 mg BSA/ml resin
Pressure flow characteristic	N/A	2 bar at 900 cm/h, 25 mm diameter \times 20 cm bed height
Maximum flow rate	5 bar at 600 cm/h, 20 cm bed height	N/A
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 1 M NaOH, 30% isopropanol or 70% ethanol. Should not be stored at < pH 3 for prolonged time	
pH stability	2 to 13	2 to 13
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 2.5 minutes residence time in 50 mM Tris-HCl, 50 mM NaCl, pH 8.0.



	BabyBio Q	OptioBio 40Q 10x100
Resin	WorkBeads 40Q	WorkBeads 40Q
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Ionic group (ligand)	Quarternary amine ($-\text{N}^+(\text{CH}_3)_3$)	Quarternary amine ($-\text{N}^+(\text{CH}_3)_3$)
Ion capacity	180 to 250 $\mu\text{mol Cl}^-/\text{ml resin}$	180 to 250 $\mu\text{mol Cl}^-/\text{ml resin}$
Dynamic binding capacity ² (DBC)	50 mg BSA/ml resin	47 mg BSA/ml resin
Column volume	1 ml 5 ml	7.9 ml
Column dimension	7 \times 28 ml (1 ml) 13 \times 38 ml (5 ml)	10 \times 100 mm
Recommended flow rate ³	BabyBio 1 ml, 0.25 to 1 ml/min (37 to 150 cm/h) BabyBio 5 ml, 1.25 to 5 ml/min (56 to 225 cm/h)	2 to 4 ml/min (150 to 300 cm/h)
Maximum flow rate ⁴	BabyBio 1 ml, 5 ml/min (780 cm/h) BabyBio 5 ml, 20 ml/min (900 cm/h)	6 ml/min (450 cm/h)
Column hardware pressure limit	0.3 MPa, 3 bar, 43 psi	2.1 MPa, 21 bar, 305 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 1 M NaOH, 30% isopropanol or 70% ethanol. Should not be stored at < pH 3 for prolonged time.	
pH stability	2 to 13	2 to 13
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4-minutes residence time in 50 mM Tris-HCl, 50 mM NaCl, pH 8.0.

³ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁴ Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature. Use half of the maximum flow rate for 20% ethanol.

	WorkBeads 40 DEAE	BabyBio DEAE
Resin	NA	WorkBeads 40 DEAE
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Ionic group (ligand)	Diethylaminoethyl ($-\text{CH}_2\text{CH}_2\text{N}^+\text{H}(\text{CH}_2\text{CH}_3)_2$)	Diethylaminoethyl ($-\text{CH}_2\text{CH}_2\text{N}^+\text{H}(\text{CH}_2\text{CH}_3)_2$)
Ion capacity	110 to 160 $\mu\text{mol Cl}^-/\text{ml resin}$	110 to 160 $\mu\text{mol Cl}^-/\text{ml resin}$
Dynamic binding capacity ² (DBC)	40 mg BSA/ml resin	40 mg BSA/ml resin
Column volume	N/A	1 ml 5 ml
Column dimension	N/A	7 \times 28 mm (1 ml) 13 \times 38 mm (5 ml)
Recommended flow rate ³	150 to 300 cm/h	BabyBio 1 ml, 0.25 to 1 ml/min (37 to 150 cm/h) BabyBio 5 ml, 1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ⁴	600 cm/h (20 cm bed height, 5 bar)	BabyBio 1 ml, 5 ml/min (780 cm/h) BabyBio 5 ml, 20 ml/min (900 cm/h)
Maximum back pressure	N/A	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 1 M NaOH, 30% isopropanol or 70% ethanol. Should not be stored at < pH 3 for prolonged time	
pH stability	3 to 13 3 to 9 (recommended pH)	3 to 13 3 to 9 (recommended pH)
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 50 mM Tris-HCl, pH 8.0.

³ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁴ Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate for 20% ethanol).



Ordering information

Product name	Pack size	Article number
WorkBeads 40Q	25 ml	40 100 001
	200 ml	40 100 002
	1 L	40 100 010
	5 L	40 100 050
	10 L	40 100 060
WorkBeads 100Q	25 ml	10 100 001
	200 ml	10 100 002
	500 ml	10 100 005
	1 L	10 100 010
	5 L	10 100 050
WorkBeads 40 DEAE	25 ml	40 150 001
	200 ml	40 150 002
	1 L	40 150 010
	5 L	40 150 050
	10 L	40 150 060
OptioBio 40Q 10x100	7.9 ml × 1	55 410 011
BabyBio IEX Screening Kit ¹	1 ml × 4	45 900 001
BabyBio Peptide Purification Kit ²	1 ml × 2	45 300 102
BabyBio Q 1 ml	1 ml × 1	45 100 101
	1 ml × 2	45 100 102
	1 ml × 5	45 100 103
	1 ml × 10	45 100 104
BabyBio Q 5 ml	5 ml × 1	45 100 105
	5 ml × 2	45 100 106
	5 ml × 5	45 100 107
	5 ml × 10	45 100 108
BabyBio DEAE 1 ml	1 ml × 1	45 150 101
	1 ml × 2	45 150 102
	1 ml × 5	45 150 103
	1 ml × 10	45 150 104
BabyBio DEAE 5 ml	5 ml × 1	45 150 105
	5 ml × 2	45 150 106
	5 ml × 5	45 150 107
	5 ml × 10	45 150 108

¹ Includes one 1 ml column of each: BabyBio S, BabyBio Q, BabyBio DEAE, BabyBio TREN.

² Bundle of: BabyBio S 1 ml × 1 and BabyBio Q 1 ml × 1.

More information

Data Sheet, DS 40 100 010

WorkBeads 40S, WorkBeads 40Q,
WorkBeads 40 DEAE

BabyBio S, BabyBio Q, BabyBio DEAE,
BabyBio TREN, BabyBio IEX Screening Kit,
BabyBio Peptide Purification Kit


www.bio-works.com/product/iex-resin
www.bio-works.com/product/optiobio-columns

Data Sheet, DS 10 200 010

WorkBeads 100S, WorkBeads 100Q

Data Sheet, DS 55 410 010

OptioBio 40S 10x100, OptioBio 40Q 10x100





Multimodal ion exchange chromatography

Multimodal ion exchange chromatography is also referred to as mixed-mode ion exchange chromatography. It utilizes the ionic interaction in combination with hydrophobic and other types of interactions. The combined effect gives the resin unique selectivities that adds new possibilities in biomolecule separation.


Target molecules

WorkBeads 40 TREN resin has a ligand that is positively charged below approx. pH 9. This resin can be used for several different applications, especially due to its higher salt tolerant properties, e.g., for alternative IEX selectivity, for sample cleanup in monoclonal antibody (mAb) purification processes to guard the protein A column from chromatins, viruses, endotoxins and other host cell impurities, or as a polishing step in the mAb purification process.



WorkBeads 40 TREN

- Differential selectivity due to higher salt tolerance and multimodal properties
- Reduced fouling of e.g. protein A resins by chromatin, viruses, endotoxins and host cell impurity removal
- High binding capacity and purity




Structure of the ligand used in WorkBeads 40 TREN

BabyBio TREN

- Prepacked for fast and reproducible purifications
- 1 ml and 5 ml columns with WorkBeads 40 TREN
- Easy to use with a syringe or chromatography system






Applications


Comparison of prepacked BabyBio TREN and BabyBio DEAE

Columns: BabyBio TREN 1 ml
BabyBio DEAE 1 ml
Binding buffer: 50 mM Tris-HCl, pH 7.4
Elution buffer: 50 mM Tris-HCl,
1 M NaCl, pH 7.4
Sample: 2.5 ml of 0.3 mg/ml
apo-transferrin,
0.2 mg/ml α -lactalbumin,
0.6 mg/ml soybean
trypsin inhibitor in
binding buffer
Flow rate: 1 ml/min (150 cm/h)
Gradient: 0 to 40% elution buffer
in 20 CV

BabyBio TREN




BabyBio DEAE



The peaks from left to right corresponds to apo-transferrin, α -lactalbumin and soybean trypsin inhibitor. The blue line corresponds to the absorbance at 280 nm and the red line to the conductivity.


Example of usage of WorkBeads 40 TREN in mAb purification processes

Flow-through mode, protection of Protein A resin (guard column).



Using WorkBeads 40 TREN as a guard column before protein A

Using WorkBeads 40 TREN upstream protein A resins is an excellent tool for eliminating the extensive bioburden on the protein A resin caused by the impurities from the host cells, and thus extending the lifetime of the protein A resin. The advantage of using WorkBeads 40 TREN upstream of WorkBeads affmAb is shown below. In this experiment up to 95% of HCP and 99% of HCD have been removed from the mAb feed loaded onto the protein A resin.



Levels of impurities in CHO cell supernatant before and after WorkBeads 40 TREN treatment. (A) HCP and (B) HCD in mAb sample loaded onto the protein A resin. (C) SDS-PAGE analyses of the feed, flow-through and eluted mAb, with or without WorkBeads 40 TREN (T) upstream WorkBeads affmAb (A).



Technical specifications

WorkBeads 40 TREN

Matrix	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm
Ligand	Tris(2-aminoethyl)amine (TAEA)
Dynamic binding capacity ²	50 mg BSA/ml resin
Maximum flow rate	600 cm/h (20 cm bed height, 5 bar)
Chemical stability	Compatible with all standard aqueous buffers used for protein purification. Should not be stored at low pH for prolonged time.
pH stability	2 to 13
Storage	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 50 mM Tris-HCl, 50 mM NaCl, pH 8.0. Optimal flow rate during binding is depending on the sample.

BabyBio TREN


Resin	WorkBeads 40 TREN
Matrix	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm
Ligand	Tris(2-aminoethyl)amine (TAEA)
Dynamic binding capacity ²	50 mg BSA/ml resin
Column volume	1 ml 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rate ³	
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ⁴	
BabyBio 1 ml	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffers used for protein purification. Do not keep the column at low pH for prolonged time.
pH stability	2 to 13
Storage	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Dynamic binding capacity determined at 4 minutes residence time (0.25 ml/min in 1 ml column) in 50 mM Tris-HCl, 50 mM NaCl, pH 8.0.

³ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

⁴ Aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate at 4°C), or by additives (e.g., use half of the maximum flow rate for 20% ethanol).



Ordering information

Product name	Pack size	Article number
WorkBeads 40 TREN	25 ml	40 603 001
	150 ml	40 603 003
	1 L	40 603 010
BabyBio TREN 1 ml	1 ml × 1	45 655 211
	1 ml × 2	45 655 212
	1 ml × 5	45 655 213
	1 ml × 10	45 655 214
BabyBio TREN 5 ml	5 ml × 1	45 655 215
	5 ml × 2	45 655 216
	5 ml × 5	45 655 217
	5 ml × 10	45 655 218
BabyBio IEX Screening Kit ¹	1 ml × 4	45 900 001

¹ Includes one 1 ml column of each: BabyBio S, BabyBio Q, BabyBio DEAE, BabyBio TREN.

More information

Data Sheet, DS 40 600 020

WorkBeads 40 TREN, BabyBio TREN

Data Sheet, DS 40 100 010

WorkBeads 40S, WorkBeads 40Q, WorkBeads 40 DEAE
 BabyBio S, BabyBio Q, BabyBio DEAE, BabyBio TREN,
 BabyBio IEX Screening Kit, BabyBio Peptide Purification Kit

www.bio-works.com/product/iex-resin






Size exclusion chromatography

Size exclusion chromatography (SEC) (also called gel filtration, GF) separates molecules on the basis of differences in size. As this is a non-binding technique the choice of running buffers is large and can be optimized for the target molecule. For best resolution when using SEC a quite slow flow rate and a sample volume of maximum 4% of the column volume should be used. This technique is therefore best suited for the polishing step in a purification process.

Target molecules

Polishing

Proteins, peptides, tagged proteins and nucleic acids. Three different pore sizes are available of WorkBeads SEC resins which make them suitable for a large range of target molecules of different sizes.



Schematic depicting size exclusion chromatography




- WorkBeads 40/100 SEC**
- WorkBeads 40/1000 SEC**
- WorkBeads 40/10 000 SEC**
- WorkBeads Macro SEC**
- WorkBeads 200 SEC**

- Produced using a proprietary cross-linking method that results in highly porous and physically stable matrices
- Availability in several different porosities give robust and wide separation ranges
- Alternative bead sizes for viscous samples
- Resistant to harsh cleaning agents (NaOH)

Comparison of WorkBeads SEC resins

	Average bead size, μm	Separation range, kD	Exclusion limit, kD	Separation range, D				
				10^4	10^5	10^6	10^7	10^8
WorkBeads 40/100 SEC	45	10 – 150	150	[Bar chart showing separation range from 10^4 to 10^5 D]				
WorkBeads 40/1000 SEC	45	10 – 1200	1200	[Bar chart showing separation range from 10^4 to 10^6 D]				
WorkBeads 40/10 000 SEC	45	10 – 10 000	10 000	[Bar chart showing separation range from 10^4 to 10^7 D]				
WorkBeads Macro SEC	45	10 – 30 000	30 000	[Bar chart showing separation range from 10^4 to $10^7.5$ D]				
WorkBeads 200 SEC	180	10 – 6000	6000	[Bar chart showing separation range from 10^4 to $10^6.5$ D]				




Applications


WorkBeads 40/1000 SEC, different flow rates

Resin: WorkBeads 40/1000 SEC
 Column: 16 × 950 mm, 181 ml
 Sample: 250 µl, 0.5 mg/ml thyroglobulin, 0.5 mg/ml ferritin, 0.5 mg/ml ovalbumin and 0.5 mg/ml ribonuclease A (in order of elution)
 Buffer: PBS, pH 7.2
 Flow rates: 25 cm/h (0.84 ml/min)
 50 cm/h (1.68 ml/min)
 100 cm/h (3.35 ml/min)


25 cm/h



50 cm/h




100 cm/h




Purification of animal viral vaccine

Resin: WorkBeads 40/1000 SEC
 Sample: Inactivated rabies virus
 (A) 10 ml (5.7% of CV)
 (B) 15 ml (8.5% of CV)
 Column: 16 × 880 mm, 176 ml
 Flow rate: 5 ml/min, 150 cm/h
 Buffer: PBS, pH 7.2

(A) Sample volume 5.7% of CV



(B) Sample volume 8.5% of CV





Technical specifications

	WorkBeads 40/100 SEC	WorkBeads 40/1000 SEC	WorkBeads 40/10 000 SEC
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Separation range ¹	10 to 150 kD	10 to 1200 kD	10 to 10 000 kD
Exclusion limit	150 kD	1200 kD	10 000 kD
Average particle size ² (D_{V50})	45 μm	45 μm	45 μm
Recommended flow rate ³	15 to 150 cm/h	15 to 150 cm/h	15 to 150 cm/h
Maximum flow rate ^{4,5}	600 cm/h	600 cm/h	300 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification. Should not be stored at low pH for prolonged time.		
pH stability	2 to 13	2 to 13	2 to 13
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ Globular proteins.

² The median particle size of the cumulative volume distribution.

³ The flow rate is important for the resolution and a lower flow rate often gives an increased resolution. A higher flow rate can be used during equilibration to speed up the separation.

⁴ Determined in water using a 25 × 200 mm column.

⁵ Note: Make sure that the column hardware max pressure is not exceeded.

	WorkBeads Macro SEC	WorkBeads 200 SEC
Separation range ¹	10 to 30 000 kD	10 to 6000 kD
Exclusion limit	30 000 kD	6000 kD
Matrix	Highly cross-linked agarose	Highly cross-linked agarose
Average particle size ² (D_{V50})	180 μm	180 μm
Recommended flow rate ³	15 to 150 cm/h	15 to 150 cm/h
Max flow rate ^{4,5}	300 cm/h	900 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification. Should not be stored at low pH for prolonged time.	
pH stability	2 to 13	2 to 13
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol


¹ Globular proteins.

² The median particle size of the cumulative volume distribution.

³ The flow rate is important for the resolution and a lower flow rate often gives an increased resolution. A higher flow rate can be used during equilibration to speed up the separation.

⁴ Determined in water using a 25 × 200 mm column.

⁵ Note: Make sure that the column hardware max pressure is not exceeded.



Ordering information

Product name	Pack size	Article number
WorkBeads 40/100 SEC	25 ml	40 340 001
	300 ml	40 340 003
	1 L	40 340 010
	5 L	40 340 050
WorkBeads 40/1000 SEC	25 ml	40 300 001
	300 ml	40 300 003
	1 L	40 300 010
	5 L	40 300 050
WorkBeads 40/10 000 SEC	25 ml	40 350 001
	300 ml	40 350 003
	1 L	40 350 010
	5 L	40 350 050
WorkBeads Macro SEC	25 ml	40 370 001
	300 ml	40 370 003
	1 L	40 370 010
	5 L	40 370 050
WorkBeads 200 SEC	300 ml	20 300 003
	1 L	20 300 010
	5 L	20 300 050


More information

Data Sheet, DS 40 300 010

WorkBeads 40/100 SEC, WorkBeads 40/1000 SEC, WorkBeads 40/10 000 SEC,
WorkBeads Macro SEC, WorkBeads 200 SEC

www.bio-works.com/product/sec-resin





Desalting/buffer exchange

Size exclusion chromatography run on low-porosity resins allows for group-separation of salt, buffer and other low molecular weight substances from larger biomolecules and proteins. This technique gives faster, simpler and a more effective desalting or buffer exchange compared to the traditional time consuming dialysis that may harm sensitive proteins and cause loss of proteins. Desalting is done with sample volumes up to 30% of the column volume and in minutes for lab scale volumes.

Target molecules

Proteins, large peptides ($M_r > 5000$), tagged proteins, nucleic acids and other biomolecules of similar size.



WorkBeads Dsalt

- Pre-swollen for fast and convenient handling
- Designed for rapid and efficient desalting and/or buffer exchange
- Group separation of high molecular weight substances from low molecular weight substances

BabyBio Dsalt


- Desalting/buffer exchange in minutes
- Keep the activity of your sensitive proteins
- Sample volumes from 20 μ l to 7.5 ml
- 1 ml and 5 ml prepacked columns
- Convenient scale-up by connecting columns in series
- Easy to use with a syringe or chromatography system



Applications

Desalting of 100 µl to 2000 µl sample


Column: BabyBio Dsalt 5 ml
 Buffer: 25 mM Na-phosphate, 150 mM NaCl, pH 7.0
 Sample: 2 mg/ml BSA in 20 mM Na-phosphate, 0.5 M NaCl, pH 7.0
 Flow rate: 5 ml/min



The solid lines correspond to absorbance at 280 nm and the dashed lines to the conductivity.

BabyBio 5 ml × 5


Column: BabyBio Dsalt 5 ml (5 columns connected in series)
 Total column volume: 25 ml
 Buffer: 25 mM Na-phosphate, 150 mM NaCl, pH 7.0
 Sample: 5 ml, 2 mg/ml BSA in 20 mM Na-phosphate, 0.5 M NaCl, pH 7.0
 Flow rate: 5 ml/min



The blue line corresponds to the absorbance at 280 nm and the red line to conductivity.

Desalting of larger sample volume

Column: WorkBeads Dsalt, 25 × 250 mm (123 ml)
 Buffer: PBS, pH 7.4
 Sample: 1 mg/ml BSA, 1 M NaCl in PBS, pH 7.4
 Sample volume: 30 ml (25 % of CV)
 37 ml (30 % of CV)
 43 ml (35% of CV)
 Flow rate: 25 ml/min (300 cm/h)



The solid lines correspond to absorbance at 280 nm (protein) and the dashed lines to the conductivity (salt).





Technical specifications

WorkBeads Dsalt

Target substance	Proteins, large peptides ($M_r > 5\,000$), nucleic acids and other biomolecules of similar size
Matrix	Highly cross-linked dextran
Average particle size ¹ (D_{v50})	150 μm
Typical sample volume	20 to 30% of the column volume (0.3 CV)
Typical flow rate	150 to 300 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 0.2 M NaOH, 0.2 M HCl, 1 M acetic acid, 8 M urea, 6 M guanidine HCl
pH stability	2 to 12
Storage	2 to 25°C in 20% ethanol or other suitable storage solution
Shipping solution	0.15% ProClin™ 150 in deionized water

¹ The median particle size of the cumulative volume distribution.

BabyBio Dsalt

Target substance	Proteins, large peptides ($M_r > 5\,000$), nucleic acids and other biomolecules of similar size
Matrix	Highly cross-linked dextran
Column volume	1 ml 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Typical sample volume	
BabyBio Dsalt 1 ml	20 to 300 μl
BabyBio Dsalt 5 ml	100 to 1500 μl
Recommended flow rate ¹	
BabyBio Dsalt 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio Dsalt 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ²	
BabyBio Dsalt 1 ml	5 ml/min (780 cm/h)
BabyBio Dsalt 5 ml	12 ml/min (540 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi
Chemical stability	Compatible with all standard aqueous buffer used for protein purification
pH stability	2 to 12
Storage	2 to 25°C in 20% ethanol

¹ Optimal flow is depending on the sample. During column wash, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively.
Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

² Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate for 20% ethanol).



Ordering information

Product name	Pack size	Article number
WorkBeads Dsalt	300 ml	40 360 003
	1 L	40 360 010
	5 L	40 360 050
	10 L	40 360 060
BabyBio Dsalt 1 ml	1 ml × 1	45 360 101
	1 ml × 2	45 360 102
	1 ml × 5	45 360 103
	1 ml × 10	45 360 104
	1 ml × 100	45 360 110
BabyBio Dsalt 5 ml	5 ml × 1	45 360 105
	5 ml × 2	45 360 106
	5 ml × 5	45 360 107
	5 ml × 10	45 360 108
	5 ml × 100	45 360 109


More information

Data Sheet, DS 40 360 010

WorkBeads Dsalt, BabyBio Dsalt

www.bio-works.com/product/sec-resin





Pre-activated resins

Pre-activated resin enables successful, convenient immobilization of ligands without the need for complex syntheses or special equipment. We have developed two different pre-activated resins where the bromohydrin active group reacts with thiol, amino and hydroxyl groups of the substance to be coupled. Two different resin porosities are available to facilitate optimized coupling of ligands of different sizes, or to optimize the prepared affinity resin for target molecules of different sizes.

Target molecules

To prepared customized chromatography resin by coupling substances with thiol, amino and hydroxyl groups.

WorkBeads 40/1000 ACT WorkBeads 40/10 000 ACT

- Ideal for coupling of specific customer designed resins
- Stable covalent linkage
- Suitable for coupling of ligands containing thiol, amino and hydroxyl groups
- Two different porosities for optimized results




BabyBio ACT

- 1 ml and 5 ml prepacked columns with WorkBeads 40/1000 ACT
- Convenient and reliable coupling procedure
- Easy to use with a syringe or chromatography system

Applications

Reaction scheme for coupling from top to bottom, primary amine, thiol and alcohol to bromohydrin activated resin.



Type of ligand and most suitable coupling conditions

Type of ligand	Functional group of ligand	Coupling conditions
Organic molecules, peptides	Thiol (Sulfhydryl) (-SH)	pH > 7 and higher
Organic molecules, peptides	Amines ¹ (-NH ₂ , -NH, -N)	pH > 8 and higher ²
Proteins, polypeptides	Thiol (Sulfhydryl) (-SH)	pH 7 and higher
Proteins, polypeptides	Primary amino (-NH ₂)	Carbonate buffer pH 8 and higher ³
Substance stable at high pH	Hydroxyl (-OH)	pH > 12 ⁴

¹ Substances containing primary, secondary and tertiary amines.

² Alkaline ligands used in excess may give high enough pH for the reaction to take place. Dissolve it in distilled water and let the basicity of the ligand determine the coupling pH.

³ Sufficient coupling without denaturation of sensitive polypeptides and proteins. Coupling reaction at a lower temperature is also possible.

⁴ High pH is required due to the low nucleophilicity of the hydroxyl group.



Technical specifications

	WorkBeads 40/1000 ACT	WorkBeads 40/10 000 ACT
Matrix	Rigid, highly cross-linked agarose	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Reactive groups	Bromohydrin	Bromohydrin
Exclusion limit	1200 kDa (globular proteins)	10 000 kDa (globular proteins)
Maximum flow rate ²	600 cm/h	600 cm/h
Reactive-groups content	200 $\mu\text{mol/ml}$	200 $\mu\text{mol/ml}$
Chemical stability (before coupling ³)	Buffers pH < 8.0	Buffers pH < 8.0
Chemical stability (after coupling ⁴)	Compatible with all standard aqueous buffers used for protein purification, 1 M NaOH, 30% isopropanol or 70% ethanol. Should not be stored at < pH 3 for prolonged time.	
pH stability ⁴	2 to 13 (after coupling)	2 to 13 (after coupling)
Storage ⁵	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Determined in water using a 10 × 300 mm column. Note: When doing a purification the optimal flow rate during binding is depending on the sample.

³ Avoid substances containing thiol and amino groups. Substances containing hydroxyl groups will only react if deprotonated.

The unreacted resin is generally stable in alcohols at neutral pH.

⁴ Agarose matrix and linker. Stability of the coupled substance may differ.

⁵ The choice of storage conditions for the coupled resin depends on the nature of the ligand.

BabyBio ACT

Resin	WorkBeads 40/1000 ACT
Matrix	Rigid, highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm
Exclusion limit	1200 kDa (globular proteins)
Reactive group	Bromohydrin
Reactive-groups content	200 $\mu\text{mol/ml}$ resin
Column volume	1 ml and 5 ml
Column dimension	7 × 28 mm (1 ml) 13 × 38 mm (5 ml)
Recommended flow rate ²	
BabyBio 1 ml	0.25 to 1 ml/min (37 to 150 cm/h)
BabyBio 5 ml	1.25 to 5 ml/min (56 to 225 cm/h)
Maximum flow rate ³	
BabyBio 1 ml	5 ml/min (780 cm/h)
BabyBio 5 ml	20 ml/min (900 cm/h)
Maximum back pressure	0.3 MPa, 3 bar, 43 psi
Chemical stability (before coupling ⁴)	Buffers pH < 8.0
Chemical stability (after coupling ⁵)	Compatible with all standard aqueous buffers used for protein purification, 1 M NaOH, 30% isopropanol or 70% ethanol. Should not be stored at < pH 3 for prolonged time.
pH stability ⁵	2 to 13 (after coupling)
Storage ⁶	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 ml/min and 5 ml/min can be used for 1 ml and 5 ml columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

³ Aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature (use half of the maximum flow rate at 4°C), or by additives (e.g. use half of the maximum flow rate for 20% ethanol).

⁴ Avoid substances containing thiol and amino groups. Substances containing hydroxyl groups will only react if deprotonated.

The unreacted resin is generally stable in alcohols at neutral pH.

⁵ Agarose matrix and linker. Stability of the coupled substance may vary.

⁶ The choice of storage conditions for the coupled resin depends on the nature of the ligand. Often 20% ethanol can be used as a bacteriostatic agent.



Ordering information

Product name	Pack size	Article number
WorkBeads 40/1000 ACT	50 ml	40 400 001
	300 ml	40 400 003
	1 L	40 400 010
	5 L	40 400 050
WorkBeads 40/10 000 ACT	50 ml	40 450 001
	300 ml	40 450 003
	1 L	40 450 010
	5 L	40 450 050
BabyBio ACT 1 ml	1 ml × 1	45 400 001
	1 ml × 2	45 400 002
	1 ml × 5	45 400 003
	1 ml × 10	45 400 004
BabyBio ACT 5 ml	5 ml × 1	45 400 005
	5 ml × 2	45 400 006
	5 ml × 5	45 400 007
	5 ml × 10	45 400 008

More information

Data Sheet, DS 40 400 010


WorkBeads 40/1000 ACT, WorkBeads 40/10 000 ACT, BabyBio ACT

www.bio-works.com/product/activated-resin









Flow rate

Binding capacity

Beads

Customer Inquiries

Mail: info@bio-works.com

Phone: +46 8 5626 7430

Web: www.bio-works.com

Company Mailing Address

Bio-Works Technologies AB
Virdings allé 18
754 50 Uppsala
Sweden

To purchase Bio-Works products
contact Bio-Works, your local
distributor or visit the web.

www.bio-works.com



*Bio-Works believes in sustainability and
care about the environment.
Please give me to a friend or recycle me.*