
DNTP Set 1 (dATP, dGTP, dCTP, dTTP)

Product Description

Product number:	UP968640
Description :	Ready-to-use dNTP solutions are available as tetrasodium salts for use in DNA polymerization reactions and all DNA labelling and sequencing processes.
Pack size :	100 mM each solution, 4 x 25 µmol each, 250 µl each vial
Biological Activity :	Not available
Reconstitution :	Not available
Stability :	Lyophilized cytokine is stable as shipped for at least 1 year at -20°C to -70 °C
Source :	Not available
Molecular weight :	dGTP Na ₄ *3H ₂ O MW 595,18 g/mol, dATP Na ₄ *3H ₂ O MW 634.2 g/mol dTTP Na ₄ *3H ₂ O MW 624.2 g/mol dCTP Na ₄ *3H ₂ O MW 591,2 g/mol
Specific activity :	Not available
Purity :	> 98% (SDS-PAGE, HPLC)

For invitro diagnostic use only/Pour usage diagnostique uniquement.