

PURIFICATION GLASSWARE




TEST TUBES

Test Tubes - Clear Borosilicate Glass EXP51

Round bottom without rim

Capacity	OD x width	Thickness	P/N	Quantity
1mL	6 x 50mm	0.5 - 0.6mm	B6AWB0	4 x 250u
4mL	10 x 75mm	0.5 - 0.6mm	AXETM5	4 x 250u
6mL	12 x 75mm	0.5 - 0.6mm	AXETN5	4 x 250u
10mL	13 x 100mm	0.5 - 0.6mm	AXETO5	4 x 250u
11mL	15 x 85mm	0.5 - 0.6mm	B6AW80	4 x 250u
15mL	16 x 100mm	0.5 - 0.6mm	AXETP5	4 x 250u
19mL	16 x 125mm	0.5 - 0.6mm	AXETR5	4 x 250u
23mL	16 x 150mm	0.5 - 0.6mm	AXETS5	4 x 250u
28.5mL	18 x 150mm	0.8mm	AXETT5	4 x 125u
36mL	20 x 150mm	0.8mm	B4FI60	4 x 125u
55mL	25 x 150mm	0.8mm	B4FI70	4 x 125u



TEST TUBES

Test tubes - Amber borosilicate Glass EXP51



Capacity	OD x width	P/N	Quantity
10mL	13 x 100mm	B4G080	1 000u
15mL	16 x 100mm	B4G090	1 000u

Test tubes - Pyrex® Glass



Manufactured from Pyrex® borosilicate glass for high resistance to thermal shock and chemical attack
Allow use at high temperature up to 500°

Capacity	OD x width	Thickness	P/N	Quantity
4mL	10 x 75mm	1.0mm	AS2Q00	100u
6mL	12 x 75mm	1.0mm	AS2Q10	100u
8mL	12 x 100mm	1.0mm	AS2Q20	100u
15mL	16 x 100mm	1.2mm	AS2Q30	100u
20mL	16 x 125mm	1.2mm	AS2Q40	100u
22mL	16 x 150mm	1.2mm	AS2Q50	100u
30mL	18 x 150mm	1.2mm	AS2Q80	100u
34mL	18 x 180mm	1.2mm	AS2Q70	100u
55mL	24 x 150mm	1.2mm	AS2Q90	100u
73mL	24 x 200mm	1.2mm	AS2QA0	50u

Test Tube Caps



Manufactured from polypropylene (autoclavable)

OD tubes (mm)	P/N	Qty
13mL	B6AWDO	1 000u
16mL	B6AWEO	1 000u
18mL	B6AWFO	1000u
20mL	B6AWGO	1 000u
25mL	B6AWHO	1 000u



CULTURE TUBES

Culture tubes - Borosilicate Glass EXP51

Without caps - round bottom

Capacity	OD x width	Thread	Thickness	P/N	Quantity
10mL	13 x 100mm	13 - 415	0.9mm	ARX3C0	1 000u
15mL	16 x 100mm	15 - 415	1.05mm	FK8505	1 000u
19mL	16 x 125mm	15 - 415	1.05mm	ARX3B0	1 000u
22mL	16 x 150mm	15 - 415	1.05mm	ARX3E0	1 000u
24mL	20 x 125mm	18 - 415	1.0mm	ARX3F0	500u

Closed screw caps (transparent polypropylene)



Thread	P/N	Qty
13 - 415	ARY3G0	1 000u
15 - 415	ARY3H0	1 000u
18 - 415	B4FIA0	500u



Closed screw caps (black phenolic)

With Rubber/PTFE seals

Thread	P/N	Qty
13 - 415	B6AX00	500u
15 - 415	B4WHY0	500u
18 - 415	B6AGK1	225u



Open screw caps (white polypropylene) with Silicone/PTFE seals

Silicone/PTFE seals

Thread	P/N	Qty
15 - 415	B2WJT0	1 000u



SOLVENT BOTTLES

Bottles, round shape



Capacity	Height x OD	Thread	Quantity	Clear Glass		Amber Glass
				Cap & pouring ring	Without cap	Without cap
10ml	54 x 36mm	GL25	10u	B688W0	B688V0	B68AS0
25ml	70 x 36mm	GL25	10u	B688Y0	B688X0	B68AV0
50ml	88 x 46mm	GL32	10u	B4REPO	B688Z0	BAREVO
100ml	100 x 56mm	GL45	10u	B689C0	B689A0	B68B00
150ml	105 x 62mm	GL45	10u	B689E0	B689D0	B68B20
250ml	138 x 70mm	GL45	10u	BV2762	B4REQ0	B68B50
500ml	176 x 86mm	GL45	10u	HP6192	B689G0	B68B80
750ml	203 x 95mm	GL45	10u	B4PQA0	B689I0	B68BA0
1 000ml	225 x 101mm	GL45	10u	HP6165	B689K0	HP6215
2 000ml	260 x 136mm	GL45	10u	HP6205	B689M0	B68BG0
3 500ml	295 x 160mm	GL45	1u	B4PQA1	B689P0	B4PQB0
5 000ml	330 x 181mm	GL45	1u	J16915	B689Q0	B68BK0
10 000ml	410 x 227mm	GL45	1u	JM5535	B689S0	B68BM0
15 000ml	465 x 258mm	GL45	1u	B689V0	B689U0	---
20 000ml	510 x 300mm	GL45	1u	B689X0	B689W0	---
500ml	153 x 101mm	GL80	10u	XF0965	---	---
1 000ml	223 x 101mm	GL80	10u	XF0975	---	---
2 000ml	510 x 300mm	GL80	10u	XF0985	---	---



SOLVENT BOTTLES

Bottles, square shape



Capacity	Height x OD	Thread	Quantity	P/N
100mL	109 x 50mm	GL32	10u	B68BWO
250mL	143 x 64mm	GL45	10u	HP6152
500mL	181 x 78mm	GL45	10u	B4REWO
1 000mL	222 x 94mm	GL45	10u	B4G2M0



B4REWO

Screw Caps



Thread	Quantity	Blue	Green	Yellow	Grey	Red caps
GL25	10u	B69N90	---	---	---	B69ND0
GL32	10u	B69NA0	---	---	---	B69NE0
GL45	10u	B4RF40	B69QR0	B69QQ0	B69QS0	B69NF0



B4RF40



B69QR0



B69QQ0



B69QS0

B69NF0



B69NE0

B69ND0

SCREW STORAGE VIALS - CONICAL BOTTOM

Vials

◆ Vials: Borosilicate glass - type 1 / Septa: High quality quality

Capacity	Height x OD	Thread	100u	250u	500u	1000u
2mL	15 x 32mm	13 - 425	---	ARYWEO	---	ARYWE1
4mL	15 x 45mm	13 - 425	---	ARYW00	---	ARYW01
8mL	16.75 x 60mm	15 - 425	---	ARYWQ0	---	ARYWQ1
10mL	19 x 70mm	22 - 400	1E4161	---	1E4162	---
20mL	26.5 x 67mm	24 - 400	1E4171	---	1E4172	---



Screw Caps

◆

Thread	Phenolic black caps - Closed Rubber / PTFE	PP white caps - Closed Silicone / PTFE	PP black caps - Open Silicone / PTFE	PP White caps - Open Silicone / PTFE	PP Black caps - Open Silicone / PTFE slit
13 - 425	ARZA00 / 250u	ARZA40 / 250u	ARZBTO / 250u	ARZAAO / 250u	ARZAGO / 250u
15 - 425	ARZA10 / 250u	ARZA50 / 250u	---	ARZABO / 250u	ARZAIO / 250u
22 - 400	ARZA20 / 100u	ARZA70 / 100u	---	ARZADO / 100u	---
24 - 400	ARZA31 / 100u	ARZA80 / 100u	---	ARZAE0 / 100u	---



Snap Storage Vials

◆

Capacity	Height x OD	Thread	Snap Vials - Clear Glass - type 3	Snap Caps - PE
			PN / 1000u	PN / 1000u
5mL	40 x 20mm	ND18	CG045A	BG846A
10mL	50 x 22mm	ND18	BG845A	BG846A
15mL	48 x 26mm	ND22	FI023A	1G057A
25mL	65 x 26mm	ND22	CV703A	1G057A



SCREW STORAGE VIALS

Vials & Kits (Vials + Caps)

● Vials: Borosilicate glass - type 1 / Septa: High quality

Capacity	Dimension	OD Caps	Vials (100u)		Unassembled kits (100u)			
			Clear vials	Amber vials	Clear vials + closed caps with PTFE seals	Amber vials + closed caps with PTFE seals	Clear vials + open caps + Silicone/ PTFE seals	Amber vials + open caps + Silicone/ PTFE seal
2mL	12 x 32mm	8mm	360190	296621	690235	690305	690370	690440
4mL	15 x 45mm	13mm	181859	360160	690249	690319	690380	690450
8mL	17 x 60mm	15mm	181869	667629	690259	690329	690390*	690460*
12mL	19 x 65mm	15mm	66763A	667649	690269	690339	690400*	690470*
16mL	21 x 70mm	18mm	181879	167470*	690279	690350*	690410*	690480*
20mL	23 x 85mm	20mm	181889	---	690289	---	690420*	---
20mL	28 x 57mm	24mm	CK1549	CE0839	OC0179	BVLOIO	---	---
30mL	28 x 73mm	24mm	CH5499	CK8259	BVLOFO	BVLOJO	---	---
40mL	28 x 95mm	24mm	217499	360259	690299	690369	690430	690500
60mL	28 x 140mm	24mm	CG0489	BX1619	BVLOGO	BVLOKO	1K1540	1K1550


*200u

Caps & Seals for Storage Vials




Thread	Closed caps with PTFE seals (100u)	Open (100u)	Silicone / PTFE seals (2,54mm thick) (100u)	Open caps with Silicone / PTFE seals assembled (100u)
8-425	97958A	282940	850550	FJ5920*
13-425	77382A	181890	438311	JQ7750*
15-425	58224A	181900	180230	B3AQFO*
18-400	58234A	181910	185010	119550
20-400	17756A	181920	176880	---
24 - 414	54594A	360210	360220	AYPNYO*

* Bonded : prevents septa fall out or push through



AUTOSAMPLER VIALS



UptiVials kits - EXP51

- ◆ Vials: Borosilicate glass - type 1
- ◆ Septa: High quality
- ◆ Dimensions: 12 x 32mm
- ◆ Regular
- ◆ Total volume: 2mL
- ◆ Usable volume: 1.5mL

UptiVials Kits (1000u)	Description	Clear glass without label	Clear glass with label
11 mm Crimp vials - 2mL	Butyl / PTFE crimp caps (aluminium)	788880	BF6490
11 mm Crimp vials - 2mL	Silicone / PTFE crimps caps (aluminium)	213910	BU0690
11 mm Snap vials - 2mL	Butyl / PTFE snap caps (PE - translucide)	1D4600	B3ASA0
11 mm Snap vials - 2mL	Silicone / PTFE snap caps (PE - translucide)	DY7080	B33U50
9mm Screw vials - 2mL	Butyl / PTFE "bonded"* screw caps (blue)	293982	744203
9mm Screw vials - 2mL	Silicone / PTFE "bonded"* screw caps (blue)	309575	T33851
9mm Screw vials - 2mL	Silicone / PTFE slit "bonded"* screw caps (blue)	T39382	CG8690
8mm Screw vials - 2mL	Butyl / PTFE screw caps (black)	B30SV0*	DY0980
8mm Screw vials - 2mL	Silicone / PTFE "bonded"* screw caps (black)	771540*	U73070
8mm Screw vials - 2mL	Silicone / PTFE slit screw caps (black)	AZ9875	AX0501

* Bonded : prevents septa fall out or push through

Tools for crimp vials 11 mm



	P/N
Hand crimper Upti-Grip	HO9150
Hand decapper Upti-Grip	HO9170



OTHERS RANGES AVAILABLE

Plastic Syringes



Syringe Filters



Safety Caps





Advion Interchim

scientific

How to request more information, a quotation or to place an order

Please see from the list below:



By Area

Europe Middle East Africa

info.EMEA@advion-interchim.com
quotes.EMEA@advion-interchim.com
orders.EMEA@advion-interchim.com
Phone: +33 4 70 03 88 55

North America

info.NAM@advion-interchim.com
quotes.NAM@advion-interchim.com
orders.NAM@advion-interchim.com
Phone: +1 607 266 9162

Asia Pacific

info.APAC@advion-interchim.com
quotes.APAC@advion-interchim.com
orders.APAC@advion-interchim.com
Phone: +46 703108410

South America

info.SAmerica@advion-interchim.com
quotes.SAmerica@advion-interchim.com
orders.SAmerica@advion-interchim.com
Phone: +46 703108410

Online

www.advion-interchim.com

All your orders will be recorded
& processed as quickly as possible.

You need technical assistance?
Our scientific experts are here to help.

Directly ask for your tailored offer

serviceconso@advion-interchim.com

Advion Interchim Scientific is the market-facing name of a partnership between Advion, Inc. and Interchim SAS, separate companies that have agreed to market and distribute their products together.
All registered and trademark belong to their respective owners. Non-contractual pictures.