

Photochemistry Starter Bundles

Everything You Need to Begin

Advances in photochemistry continue to arrive on a regular basis, enabling new and efficient pathways to key chemical synthesis reactions and techniques, including multiple C-C cross-coupling options, C-N cross-coupling, and otherwise difficult C-O cross-coupling.

EvoluChem's™ Photochemistry Starter Bundles include everything you need to further explore the opportunities inherent with photochemical reactions. Starter Bundles are available with any of our four photoreactors—PhotoRedox Box™, PhotoRedOx Box TC™, PhotoRedOx Box Duo™, and Lucent360™—either with or without sample reactions.



Seven Starter Bundle options available:

| Starter Bundle | Photoreactor | LED Light Sources | Vial Holders | Blue Light Protection | Sample Reactions |
|----------------|---------------------------------------|--|--|--|--------------------|
| HCK1020-XX-001 | PhotoRedOx Box™ HCK1006-01-016 | EvoluChem™ 365PF, 450PF HCK1012-XX-011 HCK1012-XX-002 | 8 x 4ml, 8 x 8ml HCK1006-01-019 HCK1006-01-020 | Skyper Eyewear, Orange Lens, Polycarbonate HCK1015-01-001 | 4 Sample Reactions |
| HCK1020-XX-002 | PhotoRedOx Box™ HCK1006-01-016 | EvoluChem™ 365PF, 450PF HCK1012-XX-011 HCK1012-XX-002 | 8 x 4ml, 8 x 8ml HCK1006-01-019 HCK1006-01-020 | Skyper Eyewear, Orange Lens, Polycarbonate HCK1015-01-001 | None |
| HCK1020-XX-003 | PhotoRedOx TC™ HCK1006-01-025 | EvoluChem™ 365PF, 450PF HCK1012-XX-011 HCK1012-XX-002 | 8 x 4ml, 8 x 8ml HCK1006-01-019 HCK1006-01-020 | Skyper Eyewear, Orange Lens, Polycarbonate HCK1015-01-001 | 4 Sample Reactions |
| HCK1020-XX-004 | PhotoRedOx TC™ HCK1006-01-025 | EvoluChem™ 365PF, 450PF HCK1012-XX-011 HCK1012-XX-002 | 8 x 4ml, 8 x 8ml HCK1006-01-019 HCK1006-01-020 | Skyper Eyewear, Orange Lens, Polycarbonate HCK1015-01-001 | None |
| HCK1020-XX-005 | PhotoRedOx Box Duo™ HCK1006-01-023 | EvoluChem™ 365PF, 450PF HCK1012-XX-011 HCK1012-XX-002 | 8 x 4ml (2) 8 x 8ml (2) HCK1006-01-019 HCK1006-01-020 | Skyper Eyewear, Orange Lens, Polycarbonate HCK1015-01-001 | 4 Sample Reactions |
| HCK1020-XX-006 | PhotoRedOx Box Duo™ HCK1006-01-023 | EvoluChem™ 365PF, 450PF HCK1012-XX-011 HCK1012-XX-002 | 8 x 4ml (2) 8 x 8ml (2) HCK1006-01-019 HCK1006-01-020 | Skyper Eyewear, Orange Lens, Polycarbonate HCK1015-01-001 | None |
| HCK1021-XX-020 | LUCENT360™ HCK1021-XX-001 | LUCENT360™ Side LED 450 nm (4) LUCENT360™ Bottom LED 450 nm (1) HCK1021-XX-015 HCK1021-XX-019 | 24 x 4 ml (5) HCK1021-01-004 | Skyper Eyewear, Orange Lens, Polycarbonate HCK1015-01-001 | None |

Connect seamlessly to your power source with plug options specific to your location XX = Plug Type: **01** = US (Type A); **02** = EU (Type C); **03** = UK (Type G)



© Copyright 2008-2022 HepatoChem Inc. All Rights Reserved